

Blocking Buffer 3 Sample Back-Protocol

Catalog #: KF17340

Sizes: 3X100 ml.
100 ml. General Blocker
100 ml. Neptune Blocker
100 ml. Synthetic Blocker

Storage: Block Buffers can be stored up to 12 months at 2°-8°C, and 1 week at room temperature.

All Blocking Buffers are conveniently preformulated in ready-to-use solutions that require no additional preparation.

Once your plate has been coated with an antibody or antigen, aspirate or dump out the coating solution, wash the plate 2-4 times, and simply pipette 300-400uL per well of block buffer and let it incubate.

A minimum incubation of 3 hours is recommended.

FOR RESEARCH USE ONLY

v1-01/08

Neuromics' reagents are for in vitro and certain non-human in vivo experimental use only and not intended for use in any human clinical investigation, diagnosis, prognosis, or treatment. We disclaim all liability in connection with the use of the information contained herein or otherwise, and all such risks are assumed by the user.

phone 866-350-1500 • fax 612-677-3976 • www.neuromics.com