

KAMIYA BIOMEDICAL COMPANY

Human sVE-cadherin ELISA

**For the quantitative detection of sVE-cadherin in cell culture supernatants,
human serum, plasma or other body fluids**

Cat. No. KT-031

For Research Use Only. Not for Use in Diagnostic Procedures.

PRODUCT INFORMATION

Human sVE-cadherin ELISA **Cat. No. KT-031**

PRODUCT

The **K-ASSAY®** Human sVE-cadherin ELISA is an enzyme-linked immunosorbent assay for quantitative detection of soluble human vascular endothelial (VE) cadherin in cell culture supernatants, human serum, plasma or other body fluids. The **K-ASSAY®** Human sVE-cadherin ELISA is for research use only. Not for use in diagnostic or therapeutic procedures.

DESCRIPTION

Cadherin-5, though a member of the family of cadherins, has been shown to be functionally and structurally distinct from classical cadherins (e.g. E-, N-, P-cadherins). Through its function and location cadherin-5 has been named VE-cadherin. It is a protein of a relative molecular mass of about 130 kDa.

VE-cadherin belongs to the adhesion molecules responsible for cellular interactions. The vascular endothelial cadherin (VE-cadherin) gene encodes a Ca^{2+} -dependent cell adhesion molecule required for the organization of interendothelial junctions. This gene is exclusively and constitutively expressed in endothelial cells. The corresponding protein, an endothelial-specific cadherin, is localized at the intercellular junctions. VE-cadherin mediates homophilic, calcium-dependent aggregation and cell-to-cell adhesion. In addition, it decreases intercellular permeability to high-molecular weight molecules and reduces cell migration rate across a wounded area. Thus, VE-cadherin may exert a relevant role in endothelial cell biology through control of the cohesion and organization of the intercellular junctions.

The opening of the VE-cadherin mediated endothelial barrier may be a relevant step during neutrophil extravasation. This means that despite the fact that VE-cadherin is a “nonclassical” cadherin by structure, it functions as a classic cadherin.

Vascular endothelial growth factor (VEGF) stimulation results in a maximal tyrosine phosphorylation of VE-cadherin. VE-cadherin is a transmembrane protein, the intracellular domain has been shown to interact with cytoplasmic proteins called catenins that transmit the adhesion signal upon this activation. So the VE-cadherin extracellular domain is enough for early steps of cell adhesion and recognition. However, interaction of VE-cadherin with the cytoskeleton, mediated through the cytoplasmatic domain, is necessary to provide strength and cohesion to the junction.

Apart from its established role in controlling the permeability of vascular endothelium, this molecule may have a similar role in perineurium, being important in the maintenance of the blood-nerve barrier. It furthermore functions to maintain the fibrin or collagen induced capillary tube architecture.

Specified cell adhesion molecules such as VE-cadherin are involved in the subsequent events of endothelial cell differentiation, apoptosis, and angiogenesis. In immunohistochemical studies, altered VE-cadherin expression has been described for several tumors such as haemangiomas, glioblastomas and Kaposi's sarcoma.

Most recently it has been shown that the initiation of endothelial apoptosis correlates with cleavage and disassembly of components of adherens junctions. The extracellular portion of these junctions is altered during apoptosis because VE-cadherin dramatically decreases on the surface of cells. An extracellular fragment of VE-cadherin can be detected. This shedding of VE-cadherin can be blocked by an inhibitor of metalloproteinases. It may be part of a concerted mechanism to disrupt structural and signaling properties of adherens junctions and may actively interrupt extracellular signals required for endothelial cell survival.






PRINCIPLE

An anti-VE-cadherin coating antibody is adsorbed onto microwells.

sVE-cadherin present in the sample or calibrator binds to antibodies adsorbed to the microwells; a biotin conjugated anti-VE-cadherin antibody is added and binds to sVE-cadherin captured by the first antibody.

Following incubation unbound biotin conjugated anti-VE-cadherin is removed during a wash step. Streptavidin-HRP is added and binds to the biotin conjugated anti-VE-cadherin. Following incubation unbound Streptavidin-HRP is removed during a wash step, and substrate solution reactive with HRP is added to the wells.

A colored product is formed in proportion to the amount of sVE-cadherin present in the sample. The reaction is terminated by addition of acid and absorbance is measured at 450 nm. A calibration curve is prepared from seven sVE-cadherin calibrator dilutions and sVE-cadherin sample concentration determined.



COMPONENTS

- 1 Microtiter Plate coated with Antibody to human VE-cadherin
- 1 vial (10 µL) Biotin-Conjugate concentrate anti-VE-cadherin antibody*
- 1 vial (200 µL) Streptavidin-HRP*
- 2 vials sVE-cadherin Calibrator, lyophilized, 20 ng/mL upon reconstitution
- 1 bottle (50 mL) Wash Buffer Concentrate, 20x (PBS with 1% Tween 20)*
- 1 vial (5 mL) Assay Buffer Concentrate, 20x (PBS with 1% Tween 20 and protein stabilizer)*
- 1 bottle (12 mL) Sample Diluent*
- 1 vial (7 mL) Substrate Solution I (tetramethyl-benzidine)
- 1 vial (7 mL) Substrate Solution II (0.02% buffered hydrogen peroxide)
- 1 vial (12 mL) Stop Solution (1 M Phosphoric acid)
- 1 vial (0.4 mL each) Blue-Dye, Green-Dye, Red-Dye
- 4 adhesive Plate Covers

* reagents contain preservative

Materials or Equipment Required but Not Provided

- 5 mL and 10 mL graduated pipettes
- 10 µL to 1,000 µL adjustable single channel micropipettes with disposable tips
- 50 µL to 300 µL adjustable multichannel micropipette with disposable tips
- Multichannel micropipette reservoir
- Beakers, flasks, cylinders necessary for preparation of reagents
- Device for delivery of Wash Buffer (multichannel wash bottle or automatic wash system)
- Microwell strip reader capable of reading at 450 nm (620 nm as optional reference wave length)
- Glass-distilled or deionized water
- Statistical calculator with program to perform linear regression analysis

PROTOCOLS

Specimen Collection

Cell culture supernatants, human serum, heparin plasma or other biological samples will be suitable for use in the assay. Remove serum from the clot or red cells, respectively, as soon as possible after clotting and separation.

Samples containing a visible precipitate must be clarified prior to use in the assay. Do not use grossly hemolyzed or lipemic specimens.

Clinical samples should be kept at 4°C and separated rapidly before storing at -20°C to avoid loss of bioactive VE-cadherin. If samples are to be run within 24 hours, they may be stored at 4°C. Avoid repeated freeze-thaw cycles.

For stability and suitability of samples refer to page 11.

PREPARATION OF REAGENTS

Wash Buffer

If crystals have formed in the Wash Buffer Concentrate, warm it gently until they have completely dissolved.

Pour entire contents (50 mL) of the Wash Buffer Concentrate into a clean 1,000 mL graduated cylinder. Bring final volume to 1,000 mL with glass-distilled or deionized water. Mix gently to avoid foaming. The pH of the final solution should be 7.4.

Transfer to a clean wash bottle and store at 2-25°C. Please note that the Wash Buffer is stable for 30 days. Wash Buffer may be prepared as needed according to the following table:

Number of Strips	Wash Buffer Concentrate (mL)	Distilled Water (mL)
1-6	25	475
7-12	50	950

Assay Buffer

Mix the contents of the bottle well. Add contents of Assay Buffer Concentrate (5.0 mL) to 95 mL distilled or deionized water and mix gently to avoid foaming. Store at 4°C. Please note that the Assay Buffer is stable for 30 days. Assay Buffer may be prepared as needed according to the following table:

Number of Strips	Assay Buffer Concentrate (mL)	Distilled Water (mL)
1-6	2.5	47.5
7-12	5.0	95.0

Preparation of Biotin-Conjugate

Add 90 µL of Assay Buffer to the vial containing the Biotin-Conjugate concentrate.

Make a 1:100 dilution with Assay Buffer in a clean plastic tube as needed according to the following table:

Store any diluted stock not needed immediately at -20°C.

Number of Strips	Biotin-Conjugate (mL)	Assay Buffer (mL)
1-6	0.03	2.97
7-12	0.06	5.94

Preparation of sVE-cadherin Calibrator

Reconstitute sVE-cadherin Calibrator by addition of distilled water. Reconstitution volume is stated on the label of the calibrator vial. Make sure the contents entirely dissolve by gentle swirling. The calibrator solution obtained is 20 ng/mL.

Preparation of Streptavidin-HRP

Make a 1:100 dilution of the concentrated Streptavidin-HRP solution as needed according to the following table:

Number of Strips	Streptavidin-HRP (mL)	Assay Buffer (mL)
1-6	0.06	6
7-12	0.12	12

TMB Substrate Solution

Using clean pipettes and containers known to be metal free, dispense an equal volume of Substrate Solution I into Substrate Solution II and swirl gently to mix. The TMB Substrate Solution may develop a yellow tinge over time. This does not seem to affect product performance. A blue color present in the TMB Substrate Solution, however, indicates that it has been contaminated and must be discarded.

The TMB Substrate Solution must be used within a few minutes after mixing. Warm to room temperature (RT) before use. Avoid direct exposure of TMB reagents to intense light and oxidizing agents during storage or incubation.

Substrate preparation according to assay size:

Number of Strips	Substrate Solution I (mL)	Substrate Solution II (mL)
1-6	3.0	3.0
7-12	6.0	6.0

Addition of color-giving reagents: Blue-Dye, Green-Dye, Red-Dye

In order to help our customers to avoid any mistakes in pipetting the **K-ASSAY®** Human sVE-cadherin ELISA, **KAMIYA BIOMEDICAL COMPANY** now offers a new tool that helps to monitor the addition of even very small volumes of a solution to the reaction well by giving distinctive colors to each step of the ELISA procedure.

This procedure is optional, does not in any way interfere with the test results, and is designed to help the customer with the performance of the test, but can also be omitted, just following the package insert.

Alternatively, the dye solutions from the stocks provided (Blue-Dye, Green-Dye, Red-Dye) can be added to the reagents according to the following guidelines:

1. Diluent:

Before sample dilution add the Blue-Dye at a dilution of 1:250 (see table below) to the appropriate diluent (1x) according to the test protocol. After addition of Blue-Dye, proceed according to the package insert.

5 mL Diluent	20 μ L Blue-Dye
12 mL Diluent	48 μ L Blue-Dye

2. Biotin-Conjugate:

Before dilution of the concentrated conjugate, add the Green-Dye at a dilution of 1:100 (see table below) to the Assay Buffer used for the final conjugate dilution. Proceed after addition of Green-Dye according to the package insert, preparation of Biotin-conjugate.

3 mL Assay Buffer	30 μ L Green-Dye
6 mL Assay Buffer	60 μ L Green-Dye
12 mL Assay Buffer	120 μ L Green-Dye

3. Streptavidin-HRP:

Before dilution of the concentrated Streptavidin-HRP; add the Red-Dye at a dilution of 1:250 (see table below) to the Assay Buffer used for the final Streptavidin-HRP dilution. Proceed after addition of Red-Dye according to the package insert, preparation of Streptavidin-HRP.

6 mL Assay Buffer	24 μ L Red-Dye
12 mL Assay Buffer	48 μ L Red-Dye

TEST PROTOCOL

- Mix all reagents thoroughly without foaming before use.
- Determine the number of microwell strips required to test the desired number of samples plus appropriate number of wells needed for running blanks and calibrators. Each sample, calibrator, blank, and optional control sample should be assayed in duplicate. Remove extra microwell strips coated with antibody to human VE-cadherin from holder and store in foil bag with the desiccant provided at 4°C sealed tightly.
- Wash the microwell strips twice with approximately 300 μ L Wash Buffer per well with thorough aspiration of microwell contents between washes. Take care not to scratch the surface of the microwells.
After the last wash, empty wells and tap microwell strips on absorbent pad or paper towel to remove excess Wash Buffer. Use the microwell strips immediately after washing or place upside down on a wet absorbent paper for not longer than 15 minutes. Do not allow wells to dry.
- Add 100 μ L of Sample Diluent in duplicate to all calibrator wells. Prepare calibrator dilutions by pipetting 100 μ L of reconstituted (refer to preparation of reagents) VE-cadherin Calibrator, in duplicate, into wells A1 and A2. Mix the contents of wells A1 and A2 by repeated aspiration and ejection, and transfer 100 μ L to well B1 and B2, respectively. Take care not to scratch the inner surface of the microwells. Continue this procedure five times, creating two rows of sVE-cadherin Calibrator dilutions ranging from 10 ng to 0.16 ng/mL. Discard 100 μ L of the contents from the last microwells (G1, G2) used.

Figure 1. Preparation of sVE-cadherin Calibrator dilutions:




Figure 2. Diagram depicting an example of the arrangement of blanks, calibrators and samples in the microwell strips:

	1	2	3	4
A	Calibrator 1 (10 ng/mL)	Calibrator 1 (10 ng/mL)	Sample 1	Sample 1
B	Calibrator 2 (5 ng/mL)	Calibrator 2 (5 ng/mL)	Sample 2	Sample 2
C	Calibrator 3 (2.5 ng/mL)	Calibrator 3 (2.5 ng/mL)	Sample 3	Sample 3
D	Calibrator 4 (1.25 ng/mL)	Calibrator 4 (1.25 ng/mL)	Sample 4	Sample 4
E	Calibrator 5 (0.63 ng/mL)	Calibrator 5 (0.63 ng/mL)	Sample 5	Sample 5
F	Calibrator 6 (0.32 ng/mL)	Calibrator 6 (0.32 ng/mL)	Sample 6	Sample 6
G	Calibrator 7 (0.16 ng/mL)	Calibrator 7 (0.16 ng/mL)	Sample 7	Sample 7
H	Blank	Blank	Sample 8	Sample 8

- e. Add 100 μ L of Sample Diluent in duplicate to the blank wells.
- f. Add 80 μ L of Sample Diluent, in duplicate, to the sample wells.
- g. Add 20 μ L of each Sample, in duplicate, to the designated wells.
- h. Prepare Biotin-Conjugate (refer to preparation of reagents).
- i. Add 50 μ L of diluted Biotin-Conjugate to all wells, including the blank wells.
- j. Cover with a Plate Cover and incubate at room temperature (RT) (18-25°C) for 2 hours on a microplate shaker set at 100 rpm.
- k. Remove Plate Cover and empty wells. Wash microwell strips 3 times according to step 'c' of the test protocol. Proceed immediately to the next step.
- l. Prepare Streptavidin-HRP (refer to preparation of reagents).
- m. Add 100 μ L of diluted Streptavidin-HRP to all wells, including the blank wells.
- n. Cover with a Plate Cover and incubate at RT for 1 hour on a microplate shaker at 100 rpm.
- o. Prepare TMB Substrate Solution a few minutes prior to use (refer to preparation of reagents).
- p. Remove Plate Cover and empty wells. Wash microwell strips 3 times according to step 'c' of the test protocol. Proceed immediately to the next step.
- q. Pipette 100 μ L of mixed TMB Substrate Solution to all wells, including the blank wells.
- r. Incubate the microwell strips at RT for about 10 minutes, if available on a microplate shaker at 100 rpm. Avoid direct exposure to intense light.

The color development on the plate should be monitored and substrate reaction stopped (see step 's' of the protocol) before positive wells are no longer properly recordable.

It is recommended to add the Stop Solution when the highest calibrator has developed a dark blue color.

Alternatively the color development can be monitored by the ELISA reader at 620 nm. The substrate reaction should be stopped as soon as an OD of 0.6-0.65 is reached.

- s. Stop the enzyme reaction by quickly pipetting 100 μ L of Stop Solution into each well, including the blank wells. It is important that the Stop Solution is spread quickly and uniformly throughout the microwells to completely inactivate the enzyme. Results must be read immediately after the Stop Solution is added or within one hour if the microwell strips are stored at 4°C in the dark.
- t. Read absorbance of each microwell on a spectrophotometer using 450 nm as the primary wave length (optionally 620 nm as the reference wave length; 610 nm to 650 nm is acceptable). Blank the plate reader according to the manufacturer's instructions by using the blank wells. Determine the absorbance of both, the samples and the sVE-cadherin calibrators.

Note: In case of incubation without shaking the obtained O.D. values may be lower than indicated below. Nevertheless the results are still valid.

CALCULATION OF RESULTS

- Calculate the average absorbance values for each set of duplicate calibrators and samples. Duplicates should be within 20 per cent of the mean.
- Create a calibration curve by plotting the mean absorbance for each calibrator concentration on the ordinate against the sVE-cadherin concentration on the abscissa. Draw a best fit curve through the points of the graph.
- To determine the concentration of circulating sVE-cadherin for each sample, first find the mean absorbance value on the ordinate and extend a horizontal line to the calibration curve. At the point of intersection, extend a vertical line to the abscissa and read the corresponding sVE-cadherin concentration.

For samples which have been diluted according to the instructions given in this manual 1:5, the concentration read from the calibration curve must be multiplied by the dilution factor (x5).


Note: Calculation of samples with an O.D. exceeding the range of the calibration curve may result in incorrect, low sVE-cadherin levels. Such samples should be re-analyzed at higher dilution rate in order to precisely quantitate the actual sVE-cadherin level.

It is suggested that each testing facility establishes a control sample of known sVE-cadherin concentration and runs this additional control with each assay. If the values obtained are not within the expected range of this control, the assay results may be invalid.

- A representative calibration curve is shown in Figure 3. This curve cannot be used to derive test results. Every laboratory must prepare a calibration curve for each group of microwell strips assayed.

EXPECTED RESULTS

Figure 3. Representative calibration curve for the **K-ASSAY®** Human sVE-cadherin ELISA. sVE-cadherin was diluted in serial two-fold steps in Sample Diluent, symbols represent the mean of three parallel titrations. Do not use this calibration curve to derive test results. A calibration curve must be run for each group of microwell strips assayed.



Typical data using the sVE-cadherin ELISA

Measuring wavelength: 450 nm

Reference wavelength: 620 nm

Calibrator	sVE-cadherin Concentration (ng/mL)	O.D. (450 nm)	O.D. Mean	C.V. (%)
1	10	1.954	1.983	2.0
	10	2.011		
2	5	1.101	1.083	1.6
	5	1.072		
3	2.5	0.582	0.568	3.7
	2.5	0.552		
4	1.25	0.329	0.315	6.8
	1.25	0.299		
5	0.63	0.205	0.204	1.4
	0.63	0.201		
6	0.32	0.143	0.135	9.5
	0.32	0.125		
7	0.16	0.106	0.103	4.8
	0.16	0.099		
Blank		0.080	0.083	
		0.085		

LIMITATIONS

- Since exact conditions may vary from assay to assay, a calibration curve must be established for every run.
- Bacterial or fungal contamination of either samples or reagents or cross-contamination between reagents may cause erroneous results.
- Disposable pipette tips, flasks or glassware are preferred, reusable glassware must be washed and thoroughly rinsed of all detergents before use.
- Improper or insufficient washing at any stage of the procedure will result in either false positive or false negative results. Completely empty wells before dispensing fresh Wash Buffer, fill with Wash Buffer as indicated for each wash cycle and do not allow wells to sit uncovered or dry for extended periods.

PERFORMANCE CHARACTERISTICS**Sensitivity**

The limit of detection of sVE-cadherin defined as the analyte concentration resulting in an absorption significantly higher than that of the dilution medium (mean plus two standard deviations) was determined to be less than 0.1 ng/mL (mean of 6 independent assays).

Reproducibility**Intra-assay**

Reproducibility within the assay was evaluated in three independent experiments. Each assay was carried out with 6 replicates of 8 serum samples containing different concentrations of sVE-cadherin. Two calibration curves were run on each plate. Data below show the mean sVE-cadherin concentration and the coefficient of variation for each sample. The overall intra-assay coefficient of variation has been calculated to be 4.1%.

Positive Sample	Experiment	sVE-cadherin Concentration (ng/mL)	Coefficient of Variation (%)
1	1	14.7	10.4
	2	14.0	11.1
	3	12.3	6.3
2	1	27.4	2.4
	2	33.0	4.5
	3	33.5	4.8
3	1	49.4	2.6
	2	51.5	0.7
	3	58.8	3.2
4	1	15.3	2.9
	2	16.2	4.1
	3	16.8	2.7
5	1	17.0	2.6
	2	18.1	2.1
	3	18.5	6.8
6	1	15.4	1.9
	2	16.9	2.7
	3	17.1	4.6
7	1	16.3	7.1
	2	16.6	0.9
	3	18.0	2.7
8	1	15.6	5.1
	2	16.4	3.7
	3	18.5	2.9

Inter-assay

Assay to assay reproducibility within one laboratory was evaluated in three independent experiments by three technicians. Each assay was carried out with 6 replicates of 8 serum samples containing different concentrations of sVE-cadherin. Two calibration curves were run on each plate. Data below show the mean sVE-cadherin concentration and the coefficient of variation calculated on 18 determinations of each sample. The overall inter-assay coefficient of variation has been calculated to be 7.2%.

Sample	sVE-cadherin Concentration (ng/mL)	Coefficient of Variation (%)
1	13.7	9.2
2	31.3	10.7
3	53.2	9.2
4	16.1	4.8
5	17.9	4.1
6	16.4	5.6
7	17.0	5.4
8	16.8	8.9

Spike Recovery

The spike recovery was evaluated by spiking four levels of sVE-cadherin into 4 pooled normal human serum. The amount of endogenous sVE-cadherin in unspiked serum was subtracted from the spike values. Recoveries ranged from 61% to 109% with an overall mean recovery of 82%.

Experiment	Recovery (%) sVE-cadherin Spike			
	20 ng	15 ng	10 ng	7.5 ng
1	66	69	74	61
2	77	83	86	66
3	88	93	90	78
4	99	93	109	75

Dilution Parallelism

Four serum samples with different levels of sVE-cadherin were assayed at four serial two-fold dilutions with 4 replicates each. In the table below the percent recovery of expected values is listed. Recoveries ranged from 100% to 122% with an overall mean recovery of 113%.

Sample	Dilution	sVE-cadherin Concentration (ng/mL)		
		Expected Value	Observed Value	% Recovery of Exp. Value
1	1:5	--	19.3	--
	1:10	9.7	11.3	116
	1:20	5.7	6.1	107
	1:40	3.1	3.7	119
2	1:5	--	17.6	--
	1:10	8.8	10.1	115
	1:20	5.01	5.2	103
	1:40	2.6	3.0	115
3	1:5	--	16.3	--
	1:10	8.2	9.9	120
	1:20	5.0	6.1	122
	1:40	3.0	3.5	117
4	1:5	--	25.9	--
	1:10	12.9	13.0	100
	1:20	6.5	6.8	104
	1:40	3.4	3.8	114

SAMPLE STABILITY

Freeze-Thaw Stability

Aliquots of serum and cell culture samples (unspiked or spiked) were stored frozen at -20°C and thawed up to 5 times, and sVE-cadherin levels determined. There was no significant loss of sVE-cadherin by freezing and thawing up to 5 cycles of freezing and thawing.

Storage Stability

Aliquots of a serum and cell culture samples (spiked or unspiked) were stored at -20°C, 4°C, RT and at 37°C, and the sVE-cadherin level determined after 24 hours. There was no significant loss of sVE-cadherin immunoreactivity during storage.

Comparison of Serum and Plasma

From several individuals, serum as well as heparin plasma obtained at the same time point were evaluated. sVE-cadherin concentrations were not significantly different and therefore these body fluids are suitable for the assay. It is nevertheless highly recommended to assure the uniformity of blood preparations.

Specificity

The interference of circulating factors of the immune systems was evaluated by spiking these proteins at physiologically relevant concentrations into a sVE-cadherin positive serum. There was no detectable cross reactivity.

STORAGE

Store Calibrator and Biotin-Conjugate at -20°C. Store other kit reagents at 4°C. Immediately after use remaining reagents should be returned to cold storage as indicated. Expiry of the kit and reagents is stated on labels.

The expiry of the kit components can only be guaranteed if the components are stored properly, and if, in case of repeated use of one component, the reagent is not contaminated by the first handling.

WARNINGS AND PRECAUTIONS

- Reagents are intended for research use only and are not for use in diagnostic or therapeutic procedures.
- Do not mix or substitute reagents with those from other lots or other sources.
- Do not use kit reagents beyond expiration date on label.
- Do not expose kit reagents to strong light during storage or incubation.
- Do not pipette by mouth.
- Do not eat or smoke in areas where kit reagents or samples are handled.
- Avoid contact of skin or mucous membranes with kit reagents or specimens.
- Rubber or disposable latex gloves should be worn while handling kit reagents or specimens.
- Reagents containing thimerosal as preservative may be toxic if ingested.
- Avoid contact of substrate solutions with oxidizing agents and metal.
- Avoid splashing or generation of aerosols.
- In order to avoid microbial contamination or cross-contamination of reagents or specimens which may invalidate the test use disposable pipette tips and/or pipettes.
- Use clean, dedicated reagent trays for dispensing the conjugate and substrate reagents.
- Exposure to acids will inactivate the conjugate.
- Glass-distilled water or deionized water must be used for reagent preparation.
- Substrate solutions must be at RT prior to use.
- Decontaminate and dispose specimens and all potentially contaminated materials as if they could contain infectious agents. The preferred method of decontamination is autoclaving for a minimum of 1 hour at 121.5°C.
- Liquid wastes not containing acid and neutralized waste may be mixed with sodium hypochlorite in volumes such that the final mixture contains 1.0% sodium hypochlorite. Allow 30 minutes for effective decontamination. Liquid waste containing acid must be neutralized prior to the addition of sodium hypochlorite.

REAGENT PREPARATION SUMMARY

A. Wash Buffer Add Wash Buffer Concentrate 20x (50 mL) to 950 mL distilled water.

B. Assay Buffer	Number of Strips	Assay Buffer Concentrate (mL)	Distilled Water (mL)
	1-6	2.5	47.5
	7-12	5.0	95.0

C. Biotin-Conjugate Predilution: Add 90 µL Assay Buffer to vial containing Biotin-Conjugate concentrate. Mix. Make further dilution according to the table.

	Number of Strips	prediluted 1:10 Biotin-Conjugate (mL)	Assay Buffer (mL)
	1-6	0.03	2.97
	7-12	0.06	5.94

D. sVE-cadherin Calibrator Add the volume of distilled water as stated on label to each vial of lyophilized sVE-cadherin Calibrator as needed.

E. Streptavidin-HRP

Number of Strips	Streptavidin-HRP (mL)	Assay Buffer (mL)
1-6	0.06	6
7-12	0.12	12

F. TMB Substrate Solution

Number of Strips	Substrate Solution I (mL)	Substrate Solution II (mL)
1-6	3.0	3.0
7-12	6.0	6.0

TEST PROTOCOL SUMMARY

- Wash microwell strips twice with Wash Buffer
- Add 100 μ L Sample Diluent, in duplicate, to all calibrator wells
- Pipette 100 μ L reconstituted sVE-cadherin Calibrator into the first wells and create calibrator dilutions ranging from 10 to 0.16 ng/mL by transferring 100 μ L from well to well. Discard 100 μ L from the last wells
- Add 100 μ L Sample Diluent, in duplicate, to the blank wells
- Add 80 μ L Sample Diluent to the sample wells
- Add 20 μ L Sample, in duplicate, to designated wells
- Prepare Biotin-Conjugate
- Add 50 μ L of diluted Biotin-Conjugate to all wells
- Cover microwell strips and incubate 2 hours at RT on microplate shaker
- Prepare Streptavidin-HRP
- Empty and wash microwell strips 3 times with Wash Buffer
- Add 100 μ L of diluted Streptavidin-HRP to all wells
- Cover microwell strips and incubate 1 hour at RT on microplate shaker
- Prepare TMB Substrate Solution few minutes prior to use
- Empty and wash microwell strips 3 times with Wash Buffer
- Add 100 μ L of mixed TMB Substrate Solution to all wells including blank wells
- Incubate the microwell strips for 10-20 minutes at RT on microplate shaker
- Add 100 μ L Stop Solution to all wells including blank wells
- Blank microwell reader and measure color intensity at 450 nm

Note: Calculation of samples with an O.D. exceeding 2.0 may result in incorrect, low sVE-cadherin levels. Such samples require further dilution with Sample Diluent in order to precisely quantitate the actual sVE-cadherin level.

FOR RESEARCH USE ONLY**KAMIYA BIOMEDICAL COMPANY**

12779 Gateway Drive, Seattle, WA 98168
 Tel: (206) 575-8068 Fax: (206) 575-8094
 Email: techsupport@kamyabiomedical.com
 www.kamyabiomedical.com