

OriGene Technologies Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES Phone: +1-888-267-4436 Fax: +1-301-340-8606 techsupport@origene.com

ALOOO8OPU-N OriGene EU

Acris Antibodies GmbH Schillerstr. 5 32052 Herford GERMANY Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info@acris-antibodies.com

Jejunum Lysate Catalog No.: AL00080PU-N **Quantity:** 0.5 mg **Concentration:** 2mg/ml Format: State: Liquid Purity: Jejunum tissue lysate was prepared by homogenization in lysis buffer (10 mM HEPES pH7.9, 1.5 mM MgCl₂, 10 mM KCl, 1 mM ethylenediaminetetraacetic acid, 10% glycerol, 1% NP-40, and a cocktail of protease inhibitors). Tissue and cell debris was removed by centrifugation. The product was boiled for 5 min in 1 x SDS sample buffer (50 mM Tris-HCl pH 6.8, 12.5% glycerol, 1% sodium dodecylsulfate, 0.01% bromophenol blue) containing 5% beta-mercaptoethanol. Western Blot - 10µg to 20µg per lane are recommended for mini gels. **Applications:** Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user. **Description:** Human jejunum total protein lysate Store at 2-8°C for continuous use or (in aliquots) at -70°C for longer. Storage: Avoid repeated freezing and thawing. Shelf life: 6 month from despatch.

For research and in vitro use only. Not for diagnostic or therapeutic work. Material Safety Datasheets are available at www.acris-antibodies.com or on request. Acris Antibodies is now part of the OriGene family. Learn more at www.origene.com



OG/20121030