

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES Phone: +1-888-267-4436 Fax: +1-301-340-8606 techsupport@origene.com

OriGene Technologies GmbH

Schillerstr. 5 32052 Herford GERMANY Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info-de@origene.com

AM00856PU-N

Monoclonal Antibody to Gram Positive Bacteria - Purified

Quantity:	1 mg
Concentration:	2.66 mg/ml (OD280nm, E0.1% = 1.3)
Background:	Bacteria cells are classified as Gram-positive if they retain a crystal violet dye during the Gram stain process. Gram-positive bacteria appear blue or violet under a microscope after the stain has been applied, whereas Gram-negative bacterial look red or pink. This difference in color is mainly due to the characteristics of the cell wall. Gram-positive bacteria generally have a thicker layer of peptidoglycan, a polymer consisting of sugars and amino acids that forms a homogeneous layer outside the plasma membrane. Gram-positive bacteria also have two rings supporting any flagellum and teichoic acids in the cell wall that function as as chelating agents and aid in adherence. Major groups of Gram-positive bacteria include the genera Bacillus, Listeria, Staphylococcus, Streptococcus, Enterococcus and Clostridium, as well as the phylum Actinobacteria. Gram-positive bacteria markers comprise a variety of proteins present on Gram-positive cells, and can aid in the study of function and behavior of this type of bacteria.
Host / Isotype:	Mouse / IgG1
Recommended Isotype Controls:	AM03095PU-N
Clone:	BDI813
Format:	State: Liquid purified IgG fraction (>90% pure) Purification: Protein A Chromatography Buffer System: 0.01M PBS, pH 7.2 Preservatives: 0.09% Sodium Azide Stabilizers: None
Applications:	ELISA. Immunoflourescence. Colloidal gold conjugates. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	This antibody <i>BDI813</i> clone is reactive with lipoteichoic acid (LTA) of many gram positive bacteria. Cross reacts with Listeria monocytogenes (all serotypes), Streptococcus pneumoniae, Staphylococcus aureus, Staphylococcus epidermidis, Enterococcus faecium, Bacillus cereus, Bacillus subtilis and group B Streptococcus (weak). Does <u>not</u> react with Clostridium perfringens.
Storage:	Store undiluted at 2-8°C. DO NOT FREEZE! Shelf life: one year from despatch.

For research and in vitro use only. Not for diagnostic or therapeutic work. Material Safety Datasheets are available at www.acris-antibodies.com or on request.