

OriGene Technologies Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES

Phone: +1-888-267-4436 Fax: +1-301-340-8606 techsupport@origene.com AM05532PU-T

OriGene EU

Acris Antibodies GmbH

Schillerstr. 5 32052 Herford GERMANY

Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info@acris-antibodies.com

Monoclonal Antibody to Cytokeratin 5+6+18

Alternate names: CK18, CK5, CK6, Cytokeratin 18, Cytokeratin 6, K18, K5, K6, KRT-5, KRT-6, KRT-8, KRT5, KRT6,

KRT8, Keratin type I cytoskeletal 18, Keratin type II cytoskeletal 5, Keratin type II

cytoskeletal 6

Catalog No.: AM05532PU-T

Quantity: 0.1 ml

Host / Isotype: Mouse / IgG1

Clone: LP34

Immunogen: Detergent-insoluble fraction of psoriatic human epidermis

Format: State: Liquid Tissue Culture Supernatant containing 15mM Sodium Azide

Applications: Immunohistochemistry on frozen sections: 1/100; recognises cytokeratin 5, 6 and 18 with

equal intensity, but staining of cytokeratin 18 in paraffin-embedded tissues may be

variable.

Immunohistochemistry on paraffin sections: 1/100; This product requires protein digestion

pre-treatment of paraffin sections e.g. trypsin.

Other applications not tested. Optimal dilutions are dependent on conditions and should

be determined by the user.

Specificity: This antibody recognises cytokeratins 5, 6 and 18. The antibody has a broad reactivity with

epithelial tissues ranging from simple glandular epithelial to stratified squamous epithelial. This reagent is valuable for the identification of epithelial tumours.

Species: Human.

Other species not tested.

Storage: Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

Avoid repeated freezing and thawing. Shelf life: one year from despatch.

Caution: (A full Health and Safety assessment is available upon request) This product contains

sodium azide: a POISONOUS AND HAZARDOUS SUBSTANCE which should be handled by

trained staff only.

General Readings: 1. Ghosh AK, Erber WN, Hatton CS, O'Connor NT, Falini B, Osborn M, et al. Detection of

metastatic tumour cells in routine bone marrow smears by immuno-alkaline phosphatase labelling with monoclonal antibodies. Br J Haematol. 1985 Sep;61(1):21-30. PubMed PMID:

2413878.

