

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES Phone: +1-888-267-4436 Fax: +1-301-340-8606

techsupport@origene.com

OriGene Technologies GmbH

Schillerstr. 5 32052 Herford GERMANY Phone: +49-5221-34606-0 Fax: +49-5221-34606-11

info-de@origene.com

## AM33270PU-T Monoclonal Antibody to pan Cytokeratin - Purified

Alternate names: Cytokeratin pan-reactive, pan Keratin

Quantity: 20 μg

Concentration: 0.2 mg/ml

**Background:** Cytokeratins are intermediate filament keratins found in the intracytoplasmic

cytoskeleton of epithelial tissue There are two types of Cytokeratins: the low weight, acidic type I cytokeratins and the high weight, basic or neutral type II. Cytokeratins are usually found in pairs comprising a type I Cytokeratin and a type II cytokeratin. The high molecular weight cytokeratins, which are the basic or neutral cytokeratins, comprise subtypes CK1(67), CK2(65.5), CK3(64), CK4(59), CK5(58), CK6(56), CK7(54), CK8(52.5) and CK9. The low molecular weight cytokeratins, which are the acidic cytokeratins, comprise subtypes CK10 (56.5), CK12 (56), CK13 (53), CK14 (50), CK16

(48), CK17 (46), CK18 (45), CK19 (48) and CK20 (46).

Host / Isotype: Mouse / IgG1

Recommended Isotype

**Controls:** 

SM10P (for use in human samples), SM20P (for use in rat samples), AM03095PU-N

Clone: SPM115 + SPM116

Immunogen: Human epidermal keratin.

Format: State: Liquid purified IgG fraction from Bioreactor Concentrate

Purification: Protein A/G Chromatography

Buffer System: 10mM PBS

Preservatives: 0.05% Sodium Azide

Stabilizers: 0.05% BSA

Applications: Western Blot: 0.5-1 µg/ml.

Flow Cytometry:  $0.5-1 \mu g/10^6$  cells. Immunofluorescence:  $1-2 \mu g/ml$ .

Immunohistochemistry on Formalin-Fixed Paraffin Sections: 0.5-1 µg/ml for 30

minutes at RT.

Staining of formalin-fixed tissues requires boiling tissue sections in 10mM citrate

buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 minutes.

Positive Control: Skin, Adeno- or Squamous carcinomas.

Other applications not tested. Optimal dilutions are dependent on conditions and

should be determined by the user.

**Specificity:** Twenty Human keratins are resolved with two-dimensional gel electrophoresis into

acidic (pl < 5.7) and basic (pl > 6.0) subfamilies. This antibody cocktail recognizes acidic (Type I or LMW) and basic (Type II or HMW) cytokeratins, which include CK1, CK3, CK4, CK5, CK6, CK8, CK10, CK14, CK15, CK16, and CK19. Many studies have shown

the usefulness of keratins as markers in cancer research and tumor diagnosis.

This Monoclonal Antibody is a broad spectrum anti pan-Cytokeratin antibody cocktail, which differentiates epithelial tumors from non-epithelial tumors e.g. squamous vs. adenocarcinoma of the lung, liver carcinoma, breast cancer, and esophageal cancer.



## AM33270PU-T: Monoclonal Antibody to pan Cytokeratin - Purified

It has been used to characterize the source of various neoplasms and to study the distribution of cytokeratin containing cells in epithelia during normal development

and during the development of epithelial neoplasms.

This antibody stains Cytokeratins present in normal and abnormal human tissues and has shown high sensitivity in the recognition of epithelial cells and carcinomas.

**Cellular Localization**: Cytoplasmic.

Species Reactivity: Tested: Human, Monkey, Cow, Dog, Rabbit, Mouse, Rat, Chicken.

Storage: Store undiluted at 2-8°C.

Shelf life: one year from despatch.

General Readings: 1. Woodcock-Mitchell J, Eichner R, Nelson WG, Sun TT. Immunolocalization of keratin

polypeptides in human epidermis using monoclonal antibodies. J Cell Biol. 1982

Nov;95(2 Pt 1):580-8. PubMed PMID: 6183275.

2. Tseng SCG et. al. Cell 1982; 30361.

Pictures: Formalin-Fixed, Paraffin-Embedded

Breast carcinoma stained with pan Cytokeratin Antibody cocktail Cat.-No AM33270PU (Clone SPM115+SPM116).

