

## OriGene Technologies Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES

Phone: +1-858-888-7900 Fax: +1-858-888-7904 <u>US-info@acris-antibodies.com</u>

## AP05524PU-N OriGene EU

## **Acris Antibodies GmbH**

Schillerstr. 5 32052 Herford GERMANY

Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info@acris-antibodies.com

## Polyclonal Antibody to Salmonella telhasomer - Purified

Alternate names: S. telhasomer Catalog No.: AP05524PU-N

Quantity: 1 ml

Concentration: 16.9 mg/ml
Host / Isotype: Sheep / IgG
Immunogen: Native

Format: State: Liquid Ig fraction

Purification: Ammonium sulphate precipitation

**Buffer System:** Phosphate buffered saline containing 0.09% Sodium Azide (NaN3)

**Applications:** Haemagglutination: 1/5 - 1/30.

Other applications not tested. Optimal dilutions are dependent on conditions and should

be determined by the user.

**Specificity:** The antisera has not been absorbed, and has not been tested for reactivity with somatic

(O), flagellar (H) or outer membrane polysaccharide (K) antigens.

Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

Avoid repeated freezing and thawing. Shelf life: one year from despatch.

Caution: (A full Health and Safety assessment is available upon request) This product contains

Sodium Azide: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by

trained staff only.

