

OriGene Technologies Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES Phone: +1-858-888-7900 Fax: +1-858-888-7904 <u>US-info@acris-antibodies.com</u>

AP05525PU-N OriGene EU

Acris Antibodies GmbH Schillerstr. 5 32052 Herford GERMANY Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info@acris-antibodies.com

Polyclonal Antibody to Salmonella thomasville - Purified

S. thomasville
AP05525PU-N
1 ml
20.8 mg/ml
Sheep / IgG
Native
State: Liquid Ig fraction Purification: Ammonium sulphate precipitation Buffer System: Phosphate buffered saline containing 0.09% Sodium Azide (NaN3)
Haemagglutination: 1/5 - 1/25. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
This antisera has not been absorbed, and has not been tested for reactivity with somatic (O), flagellar (H) or outer membrane polysaccharide (K) antigens.
Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. Shelf life: one year from despatch.
(A full Health and Safety assessment is available upon request) This product contains Sodium Azide: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by trained staff only.



OG/20121030