

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES Phone: +1-888-267-4436 Fax: +1-301-340-8606 techsupport@origene.com

OriGene Technologies GmbH

32052 Herford GERMANY Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info-de@origene.com

Schillerstr. 5

AP20103BT-N Polyclonal Antibody to Alcohol Oxidase - Biotin

Alternate names: AOD1, AOX, MOX, Methanol Oxidase

Quantity: 1 ml

Concentration: 10.0 mg/ml

Background: Alcohol Oxidase (AOX) is a homooctameric flavoprotein consisting of eight identical

subunits of ~74 kD, each containing a flavin adenine dinucleotide molecule (FAD) as a prosthetic group (van der Klei et al., 1991). The protein catalyzes the oxidation of methanol to formaldehyde and hydrogen peroxide, the first step in the methanol utilization pathway of certain yeasts including *Pichia pastoris*, *Hansenula*

polymorpha, and Candida boidinii (van der Klei et al., 1991). AOX is normally localized in the matrix of single membrane-bound organelles called peroxisomes. During methanol growth, the peroxisomes also contain large amounts of dihydroxyacetone synthase, the first enzyme in the methanol assimilatory pathway, and catalase (CAT), which converts the hydrogen peroxide generated by oxidases such as AOX into water and oxygen (Veenhuis and Harder, 1991). As a result, peroxisomes, which are small and few in number in glucose-grown cells, are massively induced in methanol-grown

cells (Veenhuis and Harder, 1991).

Uniprot ID: <u>000922</u> NCBI: <u>5477</u>

Host / Isotype: Rabbit / IgG

Immunogen: Alcohol dehydrogenase isolated and purified from Candida boidinii.

Freund's complete adjuvant is used in the first step of the immunization procedure.

Format: State: Lyophilized IgG fraction.

Purification: Ammonium Sulphate Precipitation and Ion Exchange Chromatography.

Buffer System: PBS, pH 7.2 without preservatives and foreign proteins.

Label: Biotin – **Conjugation Procedure:** A proprietary technique for the binding to biotin is used, followed by several purification steps. After each step activity and specificity are tested in a variety of techniques. The conjugate is lyophilized to assure

stability and long shelf life

Molar Ratio: ~4.6

Reconstitution: Restore by adding 1.0 ml of sterile distilled water.

Applications: This product is intended for use in precipitating and non-precipitating antibody-

binding assays such as e.g., ELISA and Western blotting and Immunofluorescence or

Histochemical techniques (1/1,000-1/40,000).

Other applications not tested. Optimal dilutions are dependent on conditions and

should be determined by the user.



AP20103BT-N: Polyclonal Antibody to Alcohol Oxidase - Biotin

Specificity:

The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: Immunoelectrophoresis, Cross-Immunoelectrophoresis, single Radial Immunodiffusion (Ouchterlony), block titration, ELISA, Immunoblotting and

Enzyme Inhibition.

Cross-reactivities against enzymes of other sources may occur but have not been

determined.

Recognizes Alcohol Oxidase. Species: Candida boidinii. Other species not tested.

Storage:

Store the antibody lyophilized at 2-8°C and reconstituted at 2-8°C for one week or (in

aliquots) at -20°C for longer.

If a slight precipitation occurs upon storage, this should be removed by

centrifugation.

Shelf life: one year from despatch.