

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES Phone: +1-888-267-4436 Fax: +1-301-340-8606 techsupport@origene.com **OriGene Technologies GmbH** 

Schillerstr. 5 32052 Herford GERMANY Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info-de@origene.com

## AP21390PU-N Polyclonal Antibody to Subtilisin - Aff - Purified

Alternate names:	Subtilisin E, apr, aprA, aprE, sprE
Quantity:	1 mg
Concentration:	2.0 mg/ml
Background:	Subtilisin (serine endopeptidase) is a non-specific protease (a protein-digesting enzyme) initially obtained from Bacillus subtilis. Subtilisins belong to subtilases, a group of serine proteases that initiate the nucleophilic attack on the peptide (amide) bond through a serine residue at the active site. They are physically and chemically well-characterized enzymes. Subtilisins typically have molecular weights of about 20,000 to 45,000 dalton. They can be obtained from soil bacteria, for example, <i>Bacillus amyloliquefaciens</i> . Subtilisins are secreted in large amounts from many <i>Bacillus species</i> . Subtilisin is an extracellular alkaline serine protease, it catalyzes the hydrolysis of proteins and peptide amides.
Uniprot ID:	<u>P04189</u>
NCBI:	<u>NP_388911</u>
GenelD:	<u>939313</u>
Host / Isotype:	Rabbit / IgG
Immunogen:	Subtilisin isolated and purified from <i>Bacillus subtilis</i> . Freund's complete adjuvant is used in the first step of the immunization procedure.
Format:	<ul> <li>State: Lyophilized purified IgG fraction</li> <li>Purification: Solid Phase Affinity Chromatography</li> <li>Buffer System: PBS, pH 7.2 stabilized with Dextran without preservatives and foreign proteins</li> <li>Reconstitution: Restore by adding 0.5 ml of sterile distilled water</li> </ul>
Applications:	This product is intended for use in precipitating and non-precipitating antibody- binding assays (such as e.g., ELISA and Western blotting and Immunofluorescence or Histochemical techniques), to prepare an insoluble immuno-affinity adsorbent, for labelling with a marker of choice. <u>Recommended Dilutions</u> : Non-precipitating antibody-binding techniques: 1/100-1/2,000. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.

For research and in vitro use only. Not for diagnostic or therapeutic work. Material Safety Datasheets are available at www.acris-antibodies.com or on request.

	AP21390PU-N: Polyclonal Antibody to Subtilisin - Aff - Purified
Specificity:	Subtilisin from <i>Bacillus subtilis</i> . The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: Immunoelectrophoresis, Cross-Immunoelectrophoresis, single Radial Immunodiffusion (Ouchterlony), block titration, ELISA, Immunoblotting and Enzyme Inhibition. Cross-reactivities against enzymes of other sources may occur but have not been determined. <b>Species:</b> <i>Bacillus subtilis</i> . Other species not tested.
Storage:	Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing. Shelf life: one year from despatch.

For research and in vitro use only. Not for diagnostic or therapeutic work. Material Safety Datasheets are available at www.acris-antibodies.com or on request.