

## OriGene Technologies Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES

Phone: +1-858-888-7900 Fax: +1-858-888-7904 US-info@acris-antibodies.com

## AR10715PU-N OriGene EU

## Acris Antibodies GmbH

Schillerstr. 5 32052 Herford GERMANY

Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info@acris-antibodies.com

## Human Kinase A holoenzyme type II alpha - Purified

Catalog No.: AR10715PU-N

Quantity: 3 μg

**Background:** The protein kinase A holoenzyme is a heterotetramer composed of two types of subunits:

Catalytic and Regulatory. The Catalytic subunit contains the enzyme's active site. It also contains a domain that binds ATP and a domain that binds the regulatory subunit. The Regulatory subunit consists of two molecules which bind one another in an anti-parallel orientation to form a homodimer; for type I subunits- this binding is covalent via disulfide

bonds.

This subunit also has has two domains that bind cyclic AMP, a domain that interacts with a

catalytic subunit, and an "auto-inhibitory" domain that serves as a substrate or

pseudosubstrate for the catalytic subunit.

Regulatory subunits may also have biologic activity distinct from their role in modulating catalytic subunit activity. Regulatory subunits exist in two major forms, RI and RII, with each form having two subtypes designated alpha and beta. Each of the four isotypes of the

regulatory subunit is encoded by a different gene.

In addition, three isotypes of the catalytic subunit have been identified (alpha, beta and gamma). The different isotypes tend to have different distributions within cells and among tissues. Type I enzymes inhabit cytoplasmic, soluble fractions of the cell, whereas type II

enzymes tend to associate with cellular membranes.

Species: Human
Source: E. coli

Format: State: Liquid sterile filtered clear solution in 50% Glycerol

Purity: > 95.0% as determined by SDS-PAGE

Applications: The product is suitable for analysing PKA type II agonists (cAMP analogs) or antagonists.

Other applications not tested. Optimal dilutions are dependent on conditions and should

be determined by the user.

**Description:** Inactive Holoenzyme consisting of one dimeric regulatory subunit type II alpha and two

monomeric catalytic subunits (cAMP-free).

Protein Kinase A is purified by proprietary chromatographic techniques.

Biological Activity: Holoenzyme can be activated by adding the second messenger cAMP

(Activation constant about 100nM) releasing two monomeric catalytic subunits.

Storage: Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

Avoid repeated freezing and thawing. Shelf life: one year from despatch.

