

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES Phone: +1-888-267-4436 Fax: +1-301-340-8606 techsupport@origene.com

## **OriGene Technologies GmbH**

Schillerstr. 5 32052 Herford GERMANY Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info-de@origene.com

## AR31195PU-N

## MDMA - BSA

Alternate names:	3, 4-Methylenedioxymethamphetamine, Ecstasy
Quantity:	1 mg
Concentration:	Lot Specific mg/mL (OD280nm, E 0.1% = 1.3)
Format:	<b>Purity:</b> >95% > 95% pure. Dialysis (Liquid Chromatography). <b>Buffer System:</b> 10 mM PBS, pH 7.4 <b>Preservatives:</b> 0.09% Sodium Azide
Applications:	May be used in ELISA and Lateral Flow Immunoassay. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Description:	MDMA-BSA Conjugate. 3, 4-Methylenedioxy-N-methamphetamine (MDMA) –Bovine Serum Albumin (BSA) Conjugate.
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. Shelf life: one year from despatch.

For research and in vitro use only. Not for diagnostic or therapeutic work. Material Safety Datasheets are available at www.acris-antibodies.com or on request.