

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES Phone: +1-888-267-4436 Fax: +1-301-340-8606 techsupport@origene.com

## **OriGene Technologies GmbH**

Schillerstr. 5 32052 Herford GERMANY Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info-de@origene.com

## AS00010LY-N Jurkat cells

## Jurkat cells Control Lysate

Quantity:	0.1 ml
Concentration:	Lot specific
Format:	<b>State:</b> Liquid cell lysate <b>Purity:</b> Purified <b>Buffer System:</b> Laemmli electrophoresis sample buffer
Applications:	A volume of 10-25 μl is suggested for immunoblotting. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Description:	Jurkat cells Control Lysate.
Add. Information:	Samples need to be boiled for 3-5 minutes, before loading the gel.
Storage:	Store at -20°C. Avoid repeated freezing and thawing. Shelf life: one year from despatch.

For research and in vitro use only. Not for diagnostic or therapeutic work. Material Safety Datasheets are available at www.acris-antibodies.com or on request.