



OriGene Technologies Inc.

9620 Medical Center Drive, Ste 200
Rockville, MD 20850
UNITED STATES
Phone: +1-888-267-4436
Fax: +1-301-340-8606
techsupport@origene.com

BIN129

OriGene EU

Acris Antibodies GmbH
Schillerstr. 5
32052 Herford
GERMANY
Phone: +49-5221-34606-0
Fax: +49-5221-34606-11
info@acris-antibodies.com

Echovirus 2 (Cornelis Strain) - Lysate

| | |
|--------------------------|--|
| Catalog No.: | BIN129 |
| Quantity: | 1 ml |
| Source: | TBD, BGMK cell culture Strain: Cornelis |
| Format: | State: Liquid lysate Purity: Ultra-centrifugation to remove cellular components with the resulting supernatant supplied as antigen. Buffer System: It contains no buffer or preservatives. |
| Description: | Echovirus 2 Antigen/Cornelis Strain Inactivation: Heat inactivation for 2 Hrs at 60°C. Confirmed by re-inoculating BGMK cells and no growth present. |
| Add. Information: | Internal Thread Vial. Centrifuge Product Before Opening! Heat inactivation at 2 hrs at 60 degrees Celsius. Confirmed by re-inoculating BGMK cells and no growth present. |
| Storage: | Aliquot and store at -20 to -80 degrees Celsius. Avoid multiple freeze/thaw cycles. Shelf life: six months from despatch. |
| Caution: | No test guarantees a product to be non-infectious. All materials should be handled as if potentially infectious. Generally accepted laboratory practices appropriate for infectious materials should be employed when handling this product. |