

## OriGene Technologies Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES Phone: +1-858-888-7900 Fax: +1-858-888-7904 <u>US-info@acris-antibodies.com</u>

## BIN131 OriGene EU

Acris Antibodies GmbH Schillerstr. 5 32052 Herford GERMANY Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info@acris-antibodies.com

## Enterovirus 70 (J670/71 Strain) - Lysate

Alternate names:	Enterovirus70
Catalog No.:	BIN131
Quantity:	1 ml
Source:	TBD, BGMK cell culture, Strain: J670/71
Format:	<b>State:</b> Liquid lysate. <b>Purity:</b> Ultra-centrifugation to remove cellular components with the resulting supernatant supplied as antigen.
Description:	Enterovirus 70 Antigen (J670/71 Strain) <b>Inactivation:</b> UV/Psoralen inactivation. Confirmed by re-inoculating BGMK cells and no growth present.
Storage:	Store the antigen at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid multiple freeze/thaw cycles. Shelf life: six months from despatch.
Caution:	No test guarantees a product to be non-infectious. All materials should be handled as if potentially infectious. Generally accepted laboratory practices appropriate for infectious materials should be employed when handling this product.



OG/20121030