

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES Phone: +1-888-267-4436 Fax: +1-301-340-8606 techsupport@origene.com **OriGene Technologies GmbH**

Schillerstr. 5 32052 Herford GERMANY Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info-de@origene.com

BM5502S	Monoclonal Antibody to Cytokeratin 13 - Supernatant
Alternate names:	CK-13, CK13, K13, KRT-13, KRT13, Keratin type I cytoskeletal 13, Keratin-13
Quantity:	1 ml
Concentration:	~20 µg antibody/ml
Background:	Cytokeratin 13 is a member of the keratin gene family. The keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into cytokeratins and hair keratins. Most of the type I cytokeratins consist of acidic proteins which are arranged in pairs of heterotypic keratin chains. This type I cytokeratin is paired with keratin 4 and expressed in the suprabasal layers of non- cornified stratified epithelia. Mutations in this gene and keratin 4 have been associated with the autosomal dominant disorder White Sponge Nevus.
Uniprot ID:	<u>P13646</u>
NCBI:	<u>NP_002265.2</u>
GenelD:	3860
Host / Isotype:	Mouse / IgG2b
Immunogen:	Cytokeratin purified from Human esophagus
Format:	State: Liquid Culture Supernatant Preservatives: 10mM Sodium Azide, FCS
Applications:	Immunoblotting. Immunofluorescence. Immunohistochemistry on Frozen Sections: 1/5-1/25 (±1/20), preferable in PBS. Immunohistochemistry on Paraffin Sections: 1/10 after TUF treatment. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	This Monoclonal Antibody <i>Clone 2</i> D7 does stain non-cornifying squamous epithelium. In Immunoblotting experiments the antibody recognizes only Keratin 13 (54kD) and does not show any cross reactivity with other keratins. The antibody is a useful <i>Marker</i> for histopathology. Frozen sections give excellent results in immunoperoxidase and immunofluorescence tests. Cross reactivity: Vimentin : not detectable Desmin : not detectable GFAP : not detectable Neurofilament : not detectable.
Species Reactivity:	Tested: Human.
Storage:	Store the antibody undiluted at 2-8°C. Shelf life: one year from despatch.

For research and in vitro use only. Not for diagnostic or therapeutic work. Material Safety Datasheets are available at www.acris-antibodies.com or on request.

ORIGENE BM5502S: Monoclonal Antibody to Cytokeratin 13 - Supernatant

General Readings:	 Broers JL, Ramaekers FC, Rot MK, Oostendorp T, Huysmans A, van Muijen GN, et al. Cytokeratins in different types of human lung cancer as monitored by chain-specific monoclonal antibodies. Cancer Res. 1988 Jun 1;48(11):3221-9. PubMed PMID: 2452687. van Niekerk CC, Jap PH, Thomas CM, Smeets DF, Ramaekers FC, Poels LG. Marker profile of mesothelial cells versus ovarian carcinoma cells. Int J Cancer. 1989 Jun 15;43(6):1065-71. PubMed PMID: 2659540. Ramaekers F, van Niekerk C, Poels L, Schaafsma E, Huijsmans A, Robben H, et al. Use of monoclonal antibodies to keratin 7 in the differential diagnosis of adenocarcinomas. Am J Pathol. 1990 Mar;136(3):641-55. PubMed PMID: 1690512. van Niekerk CC, Jap PH, Ramaekers FC, van de Molengraft F, Poels LG. Immunohistochemical demonstration of keratin 7 in routinely fixed paraffin- embedded human tissues. J Pathol. 1991 Oct;165(2):145-52. PubMed PMID: 1720817. van Niekerk CC, Boerman OC, Ramaekers FC, Poels LG. Marker profile of different phases in the transition of normal human ovarian epithelium to ovarian carcinomas. Am J Pathol. 1991 Feb;138(2):455-63. PubMed PMID: 1992770. Smedts F, Ramaekers F, Troyanovsky S, Pruszczynski M, Link M, Lane B, et al. Keratin expression in cervical cancer. Am J Pathol. 1992 Aug;141(2):497-511. PubMed PMID: 1379783. van Muijen GN, Ruiter DJ, Franke WW, Achtstätter T, Haasnoot WH, Ponec M, et al. Cell type heterogeneity of cytokeratin expression in complex epithelia and carcinomas as demonstrated by monoclonal antibodies specific for cytokeratins nos. 4 and 13. Exp Cell Res. 1986 Jan;162(1):97-113. PubMed PMID: 2415382. Kopan R, Traska G, Fuchs E. Retinoids as important regulators of terminal

For research and in vitro use only. Not for diagnostic or therapeutic work. Material Safety Datasheets are available at www.acris-antibodies.com or on request.