



OriGene Technologies Inc.

9620 Medical Center Drive, Ste 200  
Rockville, MD 20850  
UNITED STATES  
Phone: +1-858-888-7900  
Fax: +1-858-888-7904  
[US-info@acris-antibodies.com](mailto:US-info@acris-antibodies.com)

BP391

OriGene EU

**Acris Antibodies GmbH**  
Schillerstr. 5  
32052 Herford  
GERMANY  
Phone: +49-5221-34606-0  
Fax: +49-5221-34606-11  
[info@acris-antibodies.com](mailto:info@acris-antibodies.com)

## Polyclonal Antibody to Cortisol

**Alternate names:** Hydrocortisone

**Catalog No.:** BP391

**Quantity:** 1000 Tests

**Host:** Rabbit

**Immunogen:** Cortisol-3-BSA

**Applications:** RIA: 1/1000 - 1/2500 (125 I).

Assay method (direct liquid RIA procedure) - use as a guide only:

10-25 µl sample volume + 100 µl antibody solution + 100 µl Tracer at 1 µCi. Incubate for 90 mins at RT, then add 1 ml GaR/PEG.

Bo: 65%, ED50: 40 ng/ml.

Other applications not tested. Optimal dilutions are dependent from conditions and should be determined by the user.

Cross-Reactivity: (shown with previous batches)

Cortisol 100%

Prednisolone 36%

11-deoxycortisol 10%

Corticosterone 3.2%

Cortisone 0.9%

Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.

**Storage:** Aliquot and store the antibody at -20°C. Avoid repeated freezing and thawing. Shelf life: on

**For research and in vitro use only. Not for diagnostic or therapeutic work.**

Material Safety Datasheets are available at [www.acris-antibodies.com](http://www.acris-antibodies.com) or on request.

Acris Antibodies is now part of the OriGene family. Learn more at [www.origene.com](http://www.origene.com)



OG/20121030