

Acris Antibodies, Inc.

6815 Flanders Drive, Suite 140 San Diego, CA 92121 UNITED STATES Phone: +1-858-888-7900

Phone: +1-858-888-7900 Fax: +1-858-888-7904 US-info@acris-antibodies.com DM3054

Acris Antibodies GmbH

Schillerstr. 5 32052 Herford GERMANY

Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info@acris-antibodies.com

Monoclonal Antibody to Human Cdw32

Catalog No.: DM3054
Quantity: 0.2 mg
Concentration: 0.4 mg/ml
Host / Isotype: Mouse / IgG1

Recommended Isotype Controls:

SM10P (for use in human samples), AM03095PU-N

Clone: 7.30

Immunogen: K562 and FcgRII+L cells.

Applications: Immunofluorescence: 5-20 μg/mL. Flow cytometry: 1-5 μg/106 cells.

Immunohistochemistry: $5-20~\mu g/mL$ on frozen sections using avidin-biotin system. This antibody can be used on frozen cryostat sections as well as formalin-fixed, paraffinembedded tissue sections. Prolonged fixation is buffered formalin can destroy the epitope. For consistent and reproducible results on formalin-fixed, paraffin-embedded sections, a

high temperature antigen unmasking technique is required. This antibody was quality controlled in flow cytometry.

Recommended positive control: Macrophages, Granulocytes, Eosinophils. Other applications not tested. Optimal dilutions of this antibody are dependent on conditions

and should be determined by the user.

Other applications not tested. Optimal dilutions are dependent on conditions and should

be determined by the user.

Specificity: This antibody reacts with a 40kD (39-48kD, Endo F:35 and 29kD) polymorphic

transmembranous highly glycosylated glycoprotein, also known as Fc-gRII receptor, gp40.

Storage: Store the antibody at 4°C. Do not freeze! Shelf life: one year from despatch.

Do not dilute the entire reconstituted solution at once. Withdraw aliquots as needed with a micropipette and keep concentrated stock at 4°C. Dilute according to the particular application being used. In general, the 0.05M borate pH 8.0 containing 0.15M sodium chloride, 0.02% sodium azide, is a good dilutent to use with most antibodies. Avoid diluting the entire contents of the vial at once since the diluted solution may have reduced

stability.

General Readings: 1. lerino, B., et.al., J. Immunol. 150: 1794-1803, 1993.

2. Stuart, S.G., et.al., J. Exp. Med. 166: 1668, 1987.

3. Hogg, N., Immun

Free Call: 0800-2274746 (Germany only) - www.acris-antibodies.com

TUV NORD CERT Gmeel

TUV NORD CERT Gmeel

Roris Antihodies Gmbl