

## OriGene Technologies Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES

Phone: +1-858-888-7900 Fax: +1-858-888-7904 US-info@acris-antibodies.com

## OriGene EU

## Acris Antibodies GmbH

Schillerstr. 5 32052 Herford GERMANY

Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info@acris-antibodies.com

## C32 (human amelanotic melanoma) nuclear extract lysate (denatured)

Catalog No.:L002V3Quantity:50 μgConcentration:2

Format: Purity: Lysate quality was assured with a 12% SDS-PAGE Gel stained with Commassive

brilliant blue.

**Applications:** WB

Other applications not tested. Optimal dilutions are dependent on conditions and should

be determined by the user.

Add. Information: Lysate quality was assured with a 12% SDS-PAGE Gel stained with Commassive brilliant

blue

Storage: Store vial at -20°C or -80°C