

## OriGene Technologies Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES

Phone: +1-858-888-7900 Fax: +1-858-888-7904 <u>US-info@acris-antibodies.com</u>

## OriGene EU

## **Acris Antibodies GmbH**

Schillerstr. 5 32052 Herford GERMANY

Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info@acris-antibodies.com

## THP-1 (human acute monocytic leukemia) whole cell lysate (denatured)

Catalog No.: L007V1
Quantity: 0.2 mg
Concentration: 3

Format: Purity: Lysate quality was assured with a 12% SDS-PAGE Gel stained with Commassive

brilliant blue.

**Applications:** WB

Other applications not tested. Optimal dilutions are dependent on conditions and should

be determined by the user.

Add. Information: Cell lysate was prepared by homogenization in ice-cold modified RIPA Lysis Buffer with

cocktail of protease inhibitors (Sigma). Cell debris was removed by centrifugation. Protein concentration was determined with Bradford assay (Bio-Rad protein assay, Microplate Standard assay). The lysate was adjusted to 3.75 mg/ml and mixed with 5x SDS sample buffer to become 1xSDS sample buffer (2% SDS, 300 mM beta-mercaptoethanol, 50 mM Tris-HCl pH 6.8, 10 % Glycerol, & 0.01% Bromophenol blue) at final 3 mg/ml. The lysate was

heated at 95°C for 5 min, and cooled rapidly.

Lysate quality was assured with a 12% SDS-PAGE Gel stained with Commassive brilliant

blue.

Storage: Store vial at -20°C or -80°C

