



OriGene Technologies Inc.

9620 Medical Center Drive, Ste 200
Rockville, MD 20850
UNITED STATES
Phone: +1-858-888-7900
Fax: +1-858-888-7904
US-info@acris-antibodies.com

L022V4

OriGene EU

Acris Antibodies GmbH
Schillerstr. 5
32052 Herford
GERMANY
Phone: +49-5221-34606-0
Fax: +49-5221-34606-11
info@acris-antibodies.com

U-2 OS (human osterogenic sarcoma) nuclear extract lysate (non-denatured)

Catalog No.:	L022V4
Quantity:	50 µg
Concentration:	02. Mai
Format:	Purity: Lysate quality was assured with a 12% SDS-PAGE Gel stained with Commassive brilliant blue.
Applications:	IP, WB Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Add. Information:	Lysate quality was assured with a 12% SDS-PAGE Gel stained with Commassive brilliant blue.
Storage:	Store vial at -80°C