

Acris Antibodies GmbH

Schillerstr. 5 32052 Herford

Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info@acris-antibodies.com

293 (human embryonic kidney, adenovirus type 5 transformed) whole cell lysate (denatured)

Catalog No.: L026V1 **Quantity:** 0.2 mg **Concentration:** 3

Format: Purity: Lysate quality was assured with a 12% SDS-PAGE Gel stained with Commassive

brilliant blue.

WB **Applications:**

Other applications not tested. Optimal dilutions are dependent on conditions and should

be determined by the user.

Add. Information: Cell lysate was prepared by homogenization in ice-cold modified RIPA Lysis Buffer with

> cocktail of protease inhibitors (Sigma). Cell debris was removed by centrifugation. Protein concentration was determined with Bradford assay (Bio-Rad protein assay, Microplate Standard assay). The lysate was adjusted to 3.75 mg/ml and mixed with 5x SDS sample buffer to become 1xSDS sample buffer (2% SDS, 300 mM beta-mercaptoethanol, 50 mM Tris-HCl pH 6.8, 10 % Glycerol, & 0.01% Bromophenol blue) at final 3 mg/ml. The lysate was

heated at 95°C for 5 min, and cooled rapidly.

Lysate quality was assured with a 12% SDS-PAGE Gel stained with Commassive brilliant

blue.

Store vial at -20°C or -80°C **Storage:**

Free Call: 0800-2274746 (Germany only) - www.acris-antibodies.com