



Acris Antibodies, Inc.

6815 Flanders Drive, Suite 140
San Diego, CA 92121
UNITED STATES
Phone: +1-858-888-7900
Fax: +1-858-888-7904
US-info@acris-antibodies.com

LUM-500

Acris Antibodies GmbH

Schillerstr. 5
32052 Herford
GERMANY
Phone: +49-5221-34606-0
Fax: +49-5221-34606-11
info@acris-antibodies.com

NovaLume Chemiluminescent Western Blotting Substrate

- Catalog No.:** LUM-500
- Quantity:** 500 ml
- Background:** Novus Biologicals Chemiluminescent Substrate offers excellent sensitivity at a very economical price point.
- Applications:** , Prepare the NovaLume by mixing 1 part Reagent A with 1 part Reagent B (250 ml each). Mix well and protect from light. Completely cover the surface of the membrane with this mixture and incubate 1-5 minutes protected from light. Drain off the excess, wrap the membrane in saran and expose to film or place on imager. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
- Specificity:** NovaLume Chemiluminescent Western Blotting Substrate
- Add. Information:** Prepare the NovaLume by mixing 1 part Reagent A with 1 part Reagent B. Mix well and protect from light. Completely cover the surface of the membrane with this mixture and incubate 1-5 minutes protected from light. Drain off the excess, wrap the membrane in saran and expose to film or place on imager.
- Storage:** Store at 4C. Do not freeze.

For research and in vitro use only. Not for diagnostic or therapeutic work.

Material Safety Datasheets are available at www.acris-antibodies.com or on request.

Antibody Hotline - Technical Questions - Antibody Location Service
Free Call: 0800-2274746 (Germany only) - www.acris-antibodies.com

