



OriGene Technologies Inc.

9620 Medical Center Drive, Ste 200
Rockville, MD 20850
UNITED STATES
Phone: +1-858-888-7900
Fax: +1-858-888-7904
US-info@acris-antibodies.com

LUM-555

OriGene EU

Acris Antibodies GmbH
Schillerstr. 5
32052 Herford
GERMANY
Phone: +49-5221-34606-0
Fax: +49-5221-34606-11
info@acris-antibodies.com

NovaLume-Plus Increased Sensitivity Chemiluminescent Western Blotting Substrate

Catalog No.:	LUM-555
Quantity:	500 ml
Background:	Novus Biologicals Increased Sensitivity Chemiluminescent Substrate offers enhanced sensitivity over standard chemiluminescent reagents at a very economical price point.
Applications:	, Prepare the NovaLume-Plus by mixing 1 part Reagent A with 1 part Reagent B (250 ml each). Mix well and protect from light. Completely cover the surface of the membrane with this mixture and incubate 1-5 minutes protected from light. Drain off the excess, wrap the membrane in saran and expose to film or place on imager. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	NovaLume-Plus Increased Sensitivity Chemiluminescent Western Blotting Substrate
Add. Information:	Prepare the NovaLume-Plus by mixing 1 part Reagent A with 1 part Reagent B. Mix well and protect from light. Completely cover the surface of the membrane with this mixture and incubate 1-5 minutes protected from light. Drain off the excess, wrap the membrane in saran and expose to film or place on imager.
Storage:	Store at 4C. Do not freeze.