

# PRODUCT DATA SHEET

Anti-Nitroguanosine, clone NO<sub>2</sub>G52 **Product:** 

Cat. No.: MC-338 (50 μg)

# Background:

8-Nitroguanosine is a nitrated base of DNA and RNA. It is formed by peroxynitrite, which is generated from nitric oxide and superoxide anion radical. It is known that a large amount of nitric oxide molecules and superoxide anion, generated by inflammation, causes nitration of guanosine. Since chemically modified nucleotides cause mutation during DNA replication, 8-nitroguanosine is thought to be one of the markers of DNA damage related to mutation and cancer.

# Specificity:

Because of its very high specificity, monoclonal antibody NO<sub>2</sub>G52 recognizes 8-nitroguanine and 8-nitroguanosine, but it does not cross-react with normal nucleotide bases, 8-hydroxyguanine, 8hydroxydeoxy-guanosine, 3-nitrotyrosine, xanthine or 2-nitroimidazole. NO<sub>2</sub>G52 has very affinity for 8-nitroguanine nitroguanosine, and it slightly cross-reacts with 8-bromoguanosine, 8-bromoguanine and 8chloroguanine.

# *Immunoreactivity* [IC<sub>50</sub>]:

Strongly reacts (10 umol/l): 8-NO<sub>2</sub>-guanosine, 8-NO<sub>2</sub>-guanine.

Slightly cross-reacts (>1 mmol/l): guanosine, 8-Br-guanine, 8-Cl-guanine.

No cross-reaction: guanosine, guanine, 8-OHguanine, 8-OH-deoxyguanosine, adenine, adenosine, thymine, deoxythymidine, uracil, uridine, 3-NO<sub>2</sub>-tyrosine, 2-NO<sub>2</sub>-imidazole, cytosine.

# Species Reactivity:

All species.

## Ig Isotype:

Mouse IgG1

#### Format:

1 mg/ml solution in PBS containing 0.1% ProClin as a preservative.


# Storage:


Store at -20°C.

# Applications and Suggested Dilutions:

- ELISA: Use at 1 μg/ml
- Immunohistochemistry: Use at 10 μg/ml.

The optimal dilution for a specific application should be determined by the researcher.





Competitive ELISA using an 8-nitroguanosine-BSA-coated plate.

### Limitations:

For in vitro research use only. Not for use in diagnostics or in humans.

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.