

OriGene Technologies Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES Phone: +1-858-888-7900 Fax: +1-858-888-7904 <u>US-info@acris-antibodies.com</u>

PA1192 OriGene EU

Acris Antibodies GmbH Schillerstr. 5 32052 Herford GERMANY Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info@acris-antibodies.com

Human Menopausal Gonadotropin - Purified

Catalog No.:	PA1192
Quantity:	1 mg
Concentration:	1.0 mg/ml
Species:	Human
Source:	Urine, Urine of post-menopausal women.
Format:	 State: Sterile Filtered White lyophilized (freeze-dried) powder. Purity: >98% Greater than 98.0% as determined by: (a) Analysis by RP-HPLC. (b) Anion-exchange FPLC. (c) Analysis by reducing and non-reducing SDS-PAGE Silver Stained gel. Buffer System: Lyophilized from a concentrated solution in water containing no additives. Reconstitution: It is recommended to reconstitute the lyophilized hMG in sterile 18MΩ-cm H2O not less than 100µg/ml, which can then be further diluted to other aqueous solutions.
Description:	hMG is produced from a sterile preparation of placental glucoprotein urine of post- menopausal women. The hMG is purified by proprietary chromatographic techniques. Biological Activity: hMG conforms the following activities: 100 IU/mg of FSH and 100IU/mg LH.
Add. Information:	Protein quantitation was carried out by two independent methods: 1. UV spectroscopy at 280 nm . 2. Analysis by RP-HPLC, using a standard solution of hMG as a Reference Standard.
	Free of: HbsAg, Hepatitis B surface antigen and antibodies to HIV, Hepatitis C and HCV.
Storage:	Lyophilized hMG although stable at room temperature for 3 weeks, should be stored desiccated below -18 C. Upon reconstitution hMG should be stored at 4 C between 2-7 days and for future use below -18 C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please avoid freeze-thaw cycles.
General Readings:	 Lin YH, Seow KM, Chen HJ, Huang LW, Hwang JL, Tzeng CR. Impact of estradiol patterns in clomiphene citrate/human menopausal gonadotropin/cetrorelix protocol. Gynecol Endocrinol. 2007 Jan;23(1):45-9. PubMed PMID: 17484512. Zhang HQ, Yan B, Zhao HX, Gu DY, Jia XF, Cao L, et al. Effect of traditional Chinese herbs combined with low dose human menopausal gonadotropin applied in frozen-thawed embryo transfer. Chin J Integr Med. 2006 Dec;12(4):244-9. PubMed PMID: 17361518. Lin YH, Hwang JL, Seow KM, Huang LW, Hsieh BC, Tzeng CR. Comparison of outcome of clomiphene citrate/human menopausal gonadotropin/cetrorelix protocol and buserelin long protocola randomized study. Gynecol Endocrinol. 2006 Jun;22(6):297-302. PubMed PMID: 16785154. Zech NH, Netzbandt P, Zech H. Development of severe ovarian hyperstimulation

For research and in vitro use only. Not for diagnostic or therapeutic work. Material Safety Datasheets are available at www.acris-antibodies.com or on request.

Acris Antibodies is now part of the OriGene family. Learn more at www.origene.com



OG/20130513



PA1192: Human Menopausal Gonadotropin - Purified

syndrome after inadvertent stimulation with a gonadotropin-releasing hormone agonist and human menopausal gonadotropin in a pre-existing early pregnancy. Fertil Steril. 2005 Dec;84(6):1745. PubMed PMID: 16359980.

 Keye WR, Webster B, Dickey R, Somkuti S, Crain J, Scobey MJ. Subcutaneously administered Menopur, a new highly purified human menopausal gonadotropin, causes significantly fewer injection site reactions than Repronex in subjects undergoing in vitro fertilization. Reprod Biol Endocrinol. 2005 Nov 9;3:62. PubMed PMID: 16280073.
 Ov ulation induction in the new millennium: recombinant follicle-stimulating hormone versus human menopausal gonadotropin. Gynecol Endocrinol 2005 Mar;20(3):161-9

For research and in vitro use only. Not for diagnostic or therapeutic work. Material Safety Datasheets are available at www.acris-antibodies.com or on request. Acris Antibodies is now part of the OriGene family. Learn more at www.origene.com



OG/20130513