

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES Phone: +1-888-267-4436 Fax: +1-301-340-8606 techsupport@origene.com

OriGene Technologies GmbH

Schillerstr. 5 32052 Herford GERMANY Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info-de@origene.com

R1147 Polyclonal Antibody to Fructosyl-Amino Acid Oxidase [E.coli] Serum

Quantity: 2 ml

Concentration: 90.0 mg/ml (by Refractometry)

Background: Fructosyl amino acid oxidase (FAOX or FAOD) catalyzes the oxidative deglycation of

fructosyl amino acids, yielding the corresponding amino acid, glucosone, and hydrogen peroxide. Glucose can be attached to aminoacids by a nonenzymatic reaction called "glycation", to distinguish it from the enzymatic glycosylation of

proteins.

(FAOX) can be used for the enzymatic detection of nonenzymatically glycated proteins. Glycation of protein has been implicated in the development of diabetic complications and the aging process. Glycation of blood proteins such as hemoglobin and albumin is enhanced in diabetics with high blood glucose. The amounts of these glycated proteins reflect the level of blood glucose in periods corresponding to the half life of the protein (14 to 20 days for albumin and 1 to 2 months for hemoglobin). Since the glycation of blood proteins is not affected by transient increases in blood glucose, the levels of glycated proteins are good indices for monitoring diabetes mellitus patients during therapy. HbA1c, blood levels of which are a good index of the long term control of blood glucose in patients with diabetes mellitus, is defined as an amino end of the ß subunit (valine residue) of hemoglobin that has been glycated nonenzymatically.

Host: Sheep

Immunogen: Fructosyl-Amino Acid Oxidase [E.coli].

Format: State: Lyophilized purified Ig fraction.

Purification: Prepared from monospecific antiserum by a delipidation and

defibrination.

Buffer System: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 with

0.01% sodium azide as preservative.

Reconstitution: Restore with 2.0 ml of deionized water (or equivalent).

Applications: Suitable for Immunoblotting (Western 1:500-1:5,000), ELISA, Immunoprecipitation

(1:100) and most Immunological methods requiring high titer and specificity. Recommended Diutions: This product has been assayed against 1.0 µg of Fructosyl-Amino Acid Oxidase [E.coli] in a standard sandwich ELISA using Peroxidase conjugated Affinity Purified anti-Sheep IgG [H&L] (Rabbit) and (ABTS (2,2'-azino-bis-[3-ethylbenthiazoline-6-sulfonic acid]) as a substrate for 30 minutes at room temperature. A working dilution of 1:20,000 to 1:100,000 of the reconstitution

concentration is suggested for this product.

Other applications not tested. Optimal dilutions are dependent on conditions and

should be determined by the user.



R1147: Polyclonal Antibody to Fructosyl-Amino Acid Oxidase [E.coli] - Serum

Specificity: Assay by immunoelectrophoresis resulted in a single precipitin arc against purified

and partially purified Fructosyl-Amino Acid Oxidase [E.coli].

Cross reactivity against Fructosyl-Amino Acid Oxidase from other sources is unknown.

Storage: Store vial at 2-8°C prior to restoration. For extended storage add glycerol to 50% and

then aliquot contents and freeze at -20°C or below. Centrifuge product if not

completely clear after standing at room temperature.

This antibody is stable for one month at 2-8°C as an undiluted liquid.

Dilute only prior to immediate use. Avoid repeated freezing and thawing. Shelf life: One year from despatch.