

OriGene Technologies Inc.

9620 Medical Center Drive. Ste 200 Rockville, MD 20850 **UNITED STATES**

Phone: +1-858-888-7900 Fax: +1-858-888-7904 US-info@acris-antibodies.com

OriGene EU

R1147P

Acris Antibodies GmbH

Schillerstr. 5 32052 Herford **GERMANY**

Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info@acris-antibodies.com

Polyclonal Antibody to Fructosyl-Amino Acid Oxidase [E.coli] - Ig Fraction

R1147P Catalog No.: **Quantity:** 50 mg

Concentration: 10.0 mg/ml (by UV absorbance at 280 nm)

Fructosyl amino acid oxidase (FAOX or FAOD) catalyzes the oxidative deglycation of **Background:**

fructosyl amino acids, yielding the corresponding amino acid, glucosone, and hydrogen peroxide. Glucose can be attached to aminoacids by a nonenzymatic reaction called

"glycation", to distinguish it from the enzymatic glycosylation of proteins.

(FAOX) can be used for the enzymatic detection of nonenzymatically glycated proteins.

Glycation of protein has been implicated in the development of diabetic complications and the aging process. Glycation of blood proteins such as hemoglobin and albumin is enhanced in diabetics with high blood glucose. The amounts of these glycated proteins reflect the level of blood glucose in periods corresponding to the half life of the protein (14 to 20 days for albumin and 1 to 2 months for hemoglobin). Since the glycation of blood proteins is not affected by transient increases in blood glucose, the levels of glycated proteins are good indices for monitoring diabetes mellitus patients during therapy. HbA1c, blood levels of which are a good index of the long term control of blood glucose in patients with diabetes mellitus, is defined as an amino end of the ß subunit (valine residue) of

hemoglobin that has been glycated nonenzymatically.

Host:

Fructosyl-Amino Acid Oxidase [E.coli]. Immunogen: State: Lyophilized purified Ig fraction. Format: **Purification:** Multi-step process.

Buffer System: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 with 0.01%

(w/v) Sodium Azide as preservative.

Reconstitution: Restore with 5.0 ml of deionized water (or equivalent).

Suitable for Immunoblotting (Western or Dot blot), ELISA, Immunoprecipitation, **Applications:**

> Conjugation and most immunological methods requiring high titer and specificity. Recommended Diutions: This product has been assayed against 1.0 µg of Fructosyl-Amino

> Acid Oxidase [E.coli] in a standard sandwich ELISA using Peroxidase conjugated Affinity

Purified anti-Sheep IgG [H&L] (Goat) and ABTS (2,2'-azino-

bis-[3-ethylbenthiazoline-6-sulfonic acid]) as a substrate for 30 minutes at room

temperature.

A working dilution of 1:2,000 to 1:8,000 of the reconstitution concentration is suggested

for this product.

Other applications not tested. Optimal dilutions are dependent on conditions and should

be determined by the user.



R1147P: Polyclonal Antibody to Fructosyl-Amino Acid Oxidase [E.coli] - Ig Fraction

Specificity:

This product is an IgG fraction antibody purified from monospecific antiserum by a multistep process which includes delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Sheep Serum as well as purified and partially purified Fructosyl-Amino Acid Oxidase [E.coli]. Cross reactivity against Fructosyl-Amino Acid Oxidase from other sources is unknown.

Storage:

Store vial at 2-8°C prior to restoration. For extended storage add glycerol to 50% and then aliquot contents and freeze at -20°C or below. Centrifuge product if not completely clear after standing at room temperature.

This antibody is stable for one month at 2-8°C as an undiluted liquid.

Dilute only prior to immediate use. Avoid repeated freezing and thawing. Shelf life: One year from despatch.

