

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES Phone: +1-888-267-4436 Fax: +1-301-340-8606

techsupport@origene.com

OriGene Technologies GmbH

32052 Herford GERMANY Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info-de@origene.com

Schillerstr. 5

SA020 Human Human IgG2 (kappa chain)

Alternate names: Human Immunoglobulin G2

Quantity: 1 mg

Concentration: 1.0 mg/ml
Species: Human
Source: Myeloma
Format: State: Liquid

Purity: >95%

Buffer System: Tris buffered saline, pH 7.4, containing 0.09% sodium azide, 0.1%

EACA and 0.01% Benzamidine

Applications: ELISA.

Other applications not tested. Optimal dilutions are dependent on conditions and

should be determined by the user.

Description: Purified myeloma protein prepared from human serum

Storage: Store the IgG2 Kappa Myeloma Protein undiluted at 2-8°C for one month or (in

aliquots)

at -20°C for longer.

Avoid repeated freezing and thawing. Shelf life: one year from despatch.