

## OriGene Technologies Inc.


9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES Phone: +1-888-267-4436 Fax: +1-301-340-8606 techsupport@origene.com

## SM1667PT OriGene EU

Acris Antibodies GmbH Schillerstr. 5 32052 Herford GERMANY Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info@acris-antibodies.com

	Monoclonal Antibody to Bromodeoxyuridine / BrDU - Purified
Catalog No.:	SM1667PT
Quantity:	20 µg
Concentration:	0.5 mg/ml
Background:	The immunocytochemical detection of bromodeoxyuridine (BrdU) incorporated into DNA is a powerful tool to study the cytokinetics of normal and neoplastic cells. In vitro or in vivo labeling of tumor cells with the thymidine analogue BrdU and the subsequent detection of incorporated BrdU with specific anti-BrdU monoclonal antibodies is an accurate and comprehensive method to quantitate the degree of DNA-synthesis. BrdU is incorporated into the newly synthezised DNA of S-phase cells may provide an estimate for the fraction of cells in S-phase. Also dynamic proliferative information such as the S-phase transit rate and the potential doubling time can be obtained, by means of bivariate BrdU/DNA flow cytometric analysis.
Host / Isotype:	Rat / IgG2a
Recommended Isotype Controls:	SM15P, SM15PX
Clone:	BU1/75
Format:	State: Liquid purified IgG fraction Purification: Affinity Chromatography on Protein G Buffer System: PBS Preservatives: 0.09% Sodium Azide
Applications:	<b>Flow Cytometry:</b> Use 20 $\mu$ l of 1/25-1/200 diluted antibody to label 10 <sup>6</sup> cells in 100 $\mu$ l. Membrane permeabilization is required (See protocol below). <b>Immunohistochemistry on Formalin-Fixed, Paraffin-Embedded Sections:</b> 1/25-1/200. Denaturation of the DNA is critical for successful staining of BrdU. This can be achieved by exposing cells to heat or acid. For heat-induced epitope retrieval, 10 mM citrate buffer pH 6.0 is recommended. Alternatively, a 30 min incubation in 2M HCl can be performed. The HCl must then be neutralized for 2 min with 0.1 M Na <sub>2</sub> B <sub>4</sub> O <sub>7</sub> . Pre-treatment of tissues with proteinase K should be avoided. <b>Immunocytochemistry:</b> 1/50-1/100. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	This antibdoy reacts with BrdU in sigle stranded DNA, BrdU attached to a protein carrier or free BrdU. The antibody detects nucleated cells in S-phase with have had BrdU incorporated into their DNA.

For research and in vitro use only. Not for diagnostic or therapeutic work. Material Safety Datasheets are available at www.acris-antibodies.com or on request.





	SM1667PT: Monoclonal Antibody to Bromodeoxyuridine / BrDU - Purified
	It also reacts with Chlorodeocyuridine but with reduced staining. The antibody does not cross react with Thymidine.
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. Shelf life: one year from despatch.
Product Citations:	<b>Purchased from Acris:</b> 1. Sonal, Sidhaye J, Phatak M, Banerjee S, Mulay A, Deshpande O, et al. Myosin Vb mediated plasma membrane homeostasis regulates peridermal cell size and maintains tissue homeostasis in the zebrafish epidermis. PLoS Genet. 2014 Sep 18;10(9):e1004614. doi: 10.1371/journal.pgen.1004614. eCollection 2014 Sep. PubMed PMID: 25233349.
General Readings:	<ol> <li>Vanderlaan, M. &amp; Thomas, C.B. (1985) Characterization of monoclonal antibodies to bromodeoxyuridine. Cytometry. 6: 501-505.</li> <li>Ghiringelli, F. et al. (2005) Tumor cells convert immature myeloid dendritic cells into TGF- beta-secreting cells inducing CD4+CD25+ regulatory T cell proliferation. J. Exp. Med. 202: 919-929.</li> </ol>
	3. Dolbeare, F. (1995) Bromodeoxyuridine: a diagnostic tool in biology and medicine, Part I: Historical perspectives, histochemical methods and cell kinetics. Histochem. J. 27: 339-369.
	<ul> <li>4. Das, G. et al. (2009) Cyclin D1 fine-tunes the neurogenic output of embryonic retinal progenitor cells. Neural Dev.</li> <li>4: 15. 5. Nakhai, H. et al. (2008) Conditional ablation of Notch signaling in pancreatic</li> </ul>
	<ul> <li>development. Development. 135: 2757-65.</li> <li>6. Ghai, K. et al. (2010) Notch signaling influences neuroprotective and proliferative properties of mature Müller glia. J Neurosci. 30: 3101-12.</li> </ul>
	<ul> <li>7. Amador-Arjona, A. et al. (2011) Primary cilia regulate proliferation of amplifying progenitors in adult hippocampus: implications for learning and memory. J Neurosci. 31: 9933-44.</li> </ul>
	8. Bugler, B. et al. (2010) Unscheduled expression of CDC25B in S-phase leads to replicative stress and DNA damage. Mol Cancer. 9: 29.
	9. Gonzalo-Gobernado, R. et al (2009) Mobilization of neural stem cells and generation of new neurons in 6-OHDA-lesioned rats by intracerebroventricular infusion of liver growth factor. J Histochem Cytochem. 57: 491-502
	<ol> <li>Xu, Q. et al. (2010) Sonic hedgehog signaling confers ventral telencephalic progenitors with distinct cortical interneuron fates. Neuron. 65: 328-40.</li> <li>Zhang, J. et al. (2010) A powerful transgenic tool for fate mapping and functional</li> </ol>
	<ul> <li>analysis of newly generated neurons. BMC Neurosci. 11: 158.</li> <li>12. Bonzo, J.A. et al. (2012) Suppression of Hepatocyte Proliferation by Hepatocyte Nuclear Factor 4a in Adult Mice. J Biol Chem. Jan 12. [Epub ahead of print]</li> <li>13. Knopf, F. et al. (2011) Bone regenerates via dedifferentiation of osteoblasts in the zebrafish fin. Dev Cell. 20: 713-24.</li> </ul>
	<ul> <li>14. Grotek, B. et al. (2013) Notch signaling coordinates cellular proliferation with differentiation during zebrafish fin regeneration. Development. 140: 1412-23.</li> <li>15. Kim, T.H. et al. (2011) Genetic evidence that intestinal notch functions vary regionally and operate through a common mechanism of math1 repression. J Biol Chem. Jan 31. [Epub ahead of print]</li> </ul>
	16. Kroehne, V. et al. (2011) Regeneration of the adult zebrafish brain from neurogenic radial glia-type progenitors. Development. 138: 4831-41.
	<ol> <li>17. Liu, M.T. et al. (2009) 5-HT4 receptor-mediated neuroprotection and neurogenesis in the enteric nervous system of adult mice. J Neurosci. 29: 9683-99.</li> <li>18. Lundgren, O. et al. (2011) Intestinal epithelial stem/progenitor cells are controlled by mucosal afferent nerves. PLoS One. 6: e16295.</li> <li>19. Muja, N. et al. (2011) Neural precursors exhibit distinctly different patterns of cell</li> </ol>

**For research and in vitro use only. Not for diagnostic or therapeutic work.** Material Safety Datasheets are available at www.acris-antibodies.com or on request.

Acris Antibodies is now part of the OriGene family. Learn more at www.origene.com





## SM1667PT: Monoclonal Antibody to Bromodeoxyuridine / BrDU - Purified

migration upon transplantation during either the acute or chronic phase of EAE: A serial MR imaging study. Magn Reson Med. Feb 8. [Epub ahead of print]
20. Puverel, S. et al. (2011) RanBPM is essential for mouse spermatogenesis and oogenesis. Development. 138: 2511-21.
21. Sabo, J.K. et al. (2011) Remyelination is altered by bone morphogenic protein signaling in demyelinated lesions. J Neurosci. 31: 4504-10.
22. Scheys, J.O. et al. (2011) Evidence of adrenal failure in aging dax1-deficient mice. Endocrinology. 152: 3430-9.
23. Wang, W. et al. (2013) Extracellular signal-regulated kinase 5 (ERK5) mediates prolactin-stimulated adult neurogenesis in the subventricular zone and olfactory bulb. J Biol Chem. 288: 2623-31
24. Zemans, R.L. et al. (2013) Role of ß-catenin-Regulated CCN Matricellular Proteins in Epithelial Repair After Inflammatory Lung Injury. Am J Physiol Lung Cell Mol Physiol. Jan 11. [Epub ahead of print]

25. Zimmerman, K.M. et al. (2013) Diminished origin licensing capacity specifically sensitises tumour cells to replication stress. Mol Cancer Res. Jan 30. [Epub ahead of print] 26. Maden, M. et al. (2013) Proliferation zones in the axolotl brain and regeneration of the telencephalon. Neural Dev. 8: 1.

## Protocols: Flow Cytometry Analysis

Prepare the following solutions before proceeding: Phosphate buffered saline (PBS) 2N HCl containing 0.5% Triton X-100 PBS containing 0.05% Tween-20 PBS containing 1% BSA (PBS/BSA) 10 mg/ml Propidium iodide (PI) 0.1 M Na<sub>2</sub>B<sub>4</sub>O<sub>7</sub>, pH 8.5

As a Positive Control, BrdU labelled cells maybe obtained from Phoenix Flow Systems (http://www.phnxflow.com), catalogue number ACNC12.

1. Add BrdU to the cell suspension in culture medium to a final concentration of 10  $\mu$ mol/L and incubate for 30 min in a CO<sub>2</sub> incubator at 37°C.

2. Wash cells twice with PBS/BSA by centrifuging at 500g for 10 min, decant supernatant and resuspend in a miniumum volume of PBS.

3. Add cells slowly into 5 ml of 70% ethanol at -20°C, mixing continuously (vortex preferred). Incubate on ice for 30 min.

4. Centrifuge at 500g for 10 min, decant supernatant, and resuspend cell pellet.

5. Add 2 ml of 2N HCl containing 0.5% Triton X-100 and incubate the cells for 30 min at RT (preferably on a rocking platform).

6. Centrifuge at 500g for 10 min, decant supernatant and resuspend in 3 ml of 0.1 M  $Na_2B_4O_7$ , pH 8.5.

7. Centrifuge at 500g for 10 min, decant supernatant and resuspend the cells in PBS/BSA + 0.05% Tween-20. Adjust cell concentration to  $1 \times 10^7$ /ml.

8. Aliquot 100 ul of cell suspension into required number of 12 x 75 mm tubes.

9. Incubate the cells with the BrdU antibody at the recommended dilution for 30 min at RT.10. Add 2 ml of PBS/BSA and centrifuge the cells at 1000 rpm for 5 min.

11. If a secondary antibody layer is required then decant the supernatant and incubate the cells with the secondary antibody for 30 min at RT. If no secondary antibody layer is required then proceed to step 13.

12. Wash the cells by repeating step 10.

13. Decant off the supernatant and add 1ml of PBS containing 10  $\mu$ g/ml PI (dilute the 10 mg/ml solution of PI 1/1000 in a suitable volume of PBS)

**For research and in vitro use only. Not for diagnostic or therapeutic work.** Material Safety Datasheets are available at www.acris-antibodies.com or on request.






SM1667PT: Monoclonal Antibody to Bromodeoxyuridine / BrDU - Purified

14. Analyse cells by Flow Cytometry following the manufacturers instructions. The PI should be read on the appropriate channel set to the Peak/Area and not log scale.

**Pictures:** 

Formalin-Fixed, Paraffin-Embedded Mouse colon stained with Rat anti BrdU Antibody (Clone *BU1/75*): low power






Mouse colon stained with Rat anti BrdU Antibody Cat.-No SM1667PS (Clone *BU1/75*): low power.

Formalin-Fixed, Paraffin-Embedded

HL-60 cells were pulse labeled with BrdU for 45 minutes prior to harvesting and then incubated with primary antibody Rat anti BrdU clone BU1/75 diluted 1/100 followed by Rabbit anti Rat IgG secondary antibody (Cat.-No SP1021F) FITC-conjugated, diluted 1/200



For research and in vitro use only. Not for diagnostic or therapeutic work. Material Safety Datasheets are available at www.acris-antibodies.com or on request. Acris Antibodies is now part of the OriGene family. Learn more at www.origene.com

