

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES Phone: +1-888-267-4436 Fax: +1-301-340-8606 techsupport@origene.com

OriGene Technologies GmbH

Schillerstr. 5 32052 Herford GERMANY Phone: +49-5221-34606-0 Fax: +49-5221-34606-11

info-de@origene.com

AM02142PU-S Monoclonal Antibody to Parathyroid hormone / PTH (1-34) -

Purified

Alternate names: Parathormone, Parathyrin

Quantity: 10 μg

Background: Parathyroid hormone (PTH), or Parathormone, is secreted by the parathyroid glands

as a polypeptide containing 84 amino acids. It acts to increase the concentration of calcium in the blood, whereas calcitonin (a hormone produced by the parafollicular

cells of the thyroid gland) acts to decrease calcium concentration.

Uniprot ID: P01270

NCBI: <u>NP 000306.1</u>

GeneID: <u>5741</u>

Host / Isotype: Rat / IgG2b

Clone: 3B3

Immunogen: Synthetic Human PTH (aa 1-34) BSA conjugated.

Genename: PTH AA Sequence:

SVSEIQLMHNLGKHLNSMERVEWLRKKLQDVHNF

Format: State: Lyophilized purified Cell supernatant

Purification: Protein G Chromatography

Buffer System: PBS, pH 7.4 without preservatives **Reconstitution:** Restore in aqua bidest to 1 mg/ml.

Applications: ELISA.

Immunoprecipitation (Affinity Chromatography) (Ref. 1).

Immunohistochemistry on CryoSections and Paraffin-Sections: 1/25-50.

Other applications not tested. Optimal dilutions are dependent on conditions and

should be determined by the user.

Specificity: This antibody detects Human oxidized PTH (aa 1-34), native PTH in tissue sections and

in circulation.

Cross reacts with native PTH (aa 1-84).

Species: Human.

Other species not tested.

Add. Information: LocusID 5741

Storage: Store lyophilized at 2-8°C for 6 months or at -20°C long term.

After reconstitution store the antibody undiluted at 2-8°C for one month

or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing. Shelf life: one year from despatch.

AM02142PU-S: Monoclonal Antibody to Parathyroid hormone / PTH (1-34) - Purified



General Readings:

1. Hocher B, Armbruster FP, Stoeva S, Reichetzeder C, Grön HJ, Lieker I, et al. Measuring parathyroid hormone (PTH) in patients with oxidative stress--do we need a fourth generation parathyroid hormone assay? PLoS One. 2012;7(7):e40242. doi: 10.1371/journal.pone.0040242. Epub 2012 Jul 6. PubMed PMID: 22792251.

2. Logue FC, Perry B, Biggart EM, Chapman RS, Beastall GH. Production and characterisation of monoclonal antibodies to parathyroid hormone (1-34). J Immunol Methods. 1991 Mar 21;137(2):159-66. PubMed PMID: 2013695.