

OriGene Technologies Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES Phone: +1-858-888-7900 Fax: +1-858-888-7904 <u>US-info@acris-antibodies.com</u>

AP05513PU-N OriGene EU

Acris Antibodies GmbH Schillerstr. 5 32052 Herford GERMANY Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info@acris-antibodies.com

Polyclonal Antibody to Salmonella give - Purified

| Alternate names: | S. give |
|-----------------------|---|
| Catalog No.: | AP05513PU-N |
| Quantity: | 1 ml |
| Concentration: | 29.8 mg/ml |
| Host / Isotype: | Sheep / IgG |
| Immunogen: | Native |
| Format: | State: Liquid Ig fraction Purification: Ammonium sulphate precipitation Buffer System: Phosphate buffered saline containing 0.09% Sodium Azide (NaN3) |
| Applications: | Haemagglutination: 1/10 - 1/50. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user. |
| Specificity: | Not tested for reactivity with somatic (O), flagellar (H) or outer membrane polysaccharide (K) antigens. |
| Storage: | Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. Shelf life: one year from despatch. |
| Caution: | (A full Health and Safety assessment is available upon request) This product contains Sodium Azide: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by trained staff only. |



OG/20121030