



## OriGene Technologies Inc.

9620 Medical Center Drive, Ste 200  
Rockville, MD 20850  
UNITED STATES  
Phone: +1-858-888-7900  
Fax: +1-858-888-7904  
[US-info@acris-antibodies.com](mailto:US-info@acris-antibodies.com)

AP05519PU-N

## OriGene EU

**Acris Antibodies GmbH**  
Schillerstr. 5  
32052 Herford  
GERMANY  
Phone: +49-5221-34606-0  
Fax: +49-5221-34606-11  
[info@acris-antibodies.com](mailto:info@acris-antibodies.com)

# Polyclonal Antibody to Salmonella paratyphi C - Purified

|                         |   |
|-------------------------|---|
| <b>Alternate names:</b> | S. paratyphi C  |
| <b>Catalog No.:</b>     | AP05519PU-N   |
| <b>Quantity:</b>        | 1 ml  |
| <b>Concentration:</b>   | 20.0 mg/ml  |
| <b>Host / Isotype:</b>  | Sheep / IgG   |
| <b>Immunogen:</b>       | Native  |
| <b>Format:</b>          | <b>State:</b> Liquid Ig fraction<br><b>Purification:</b> Ammonium sulphate precipitation<br><b>Buffer System:</b> Phosphate buffered saline containing 0.09% Sodium Azide (NaN <sub>3</sub> ) |
| <b>Applications:</b>    | ELISA.<br>Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.  |
| <b>Specificity:</b>     | This antibody reacts to Group C1 <i>Salmonella paratyphi C</i>  |
| <b>Storage:</b>         | Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.<br>Avoid repeated freezing and thawing.<br>Shelf life: one year from despatch.                      |
| <b>Caution:</b>         | (A full Health and Safety assessment is available upon request) This product contains Sodium Azide: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by trained staff only.       |