

OriGene Technologies Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES Phone: +1-858-888-7900 Fax: +1-858-888-7904 <u>US-info@acris-antibodies.com</u>

AP05521PU-N OriGene EU

Acris Antibodies GmbH Schillerstr. 5 32052 Herford GERMANY Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info@acris-antibodies.com

Polyclonal Antibody to Salmonella seattle - Purified

Alternate names:	S. seattle
Catalog No.:	AP05521PU-N
Quantity:	1 ml
Concentration:	17.9 mg/ml
Host / Isotype:	Sheep / IgG
Immunogen:	Native
Format:	State: Liquid Ig fraction Purification: Ammonium sulphate precipitation Buffer System: Phosphate buffered saline containing 0.09% Sodium Azide (NaN3)
Applications:	Haemagglutination: 1/5 - 1/30. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	This antisera has not been absorbed, and has not been tested for reactivity with somatic (O), flagellar (H) or outer membrane polysaccharide (K) antigens.
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing. Shelf life: one year from despatch.
Caution:	(A full Health and Safety assessment is available upon request) This product contains Sodium Azide: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by trained staff only.



OG/20121030