

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES Phone: +1-888-267-4436 Fax: +1-301-340-8606 techsupport@origene.com

OriGene Technologies GmbH

Schillerstr. 5 32052 Herford GERMANY Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info-de@origene.com

AP08139SU-N Normal Rabbit Serum

Quantity: 5 ml

Background: Serum is the clear, straw-colored, liquid portion of blood plasma that does not

contain fibrinogen or blood cells and remains fluid after clotting.

Serums are obtained from non-hemolyzed blood that is collected from healthy and fasted donors. The donors are not medicated in any way and are maintained on an antibiotic free diet. Sterilization is accomplished by Millipore filtration. These serums

provide excellent growth promoting components for tissue culture and

microbiological organisms.

Host: Rabbit

Format: State: Serum (Sterile filtered) containing 0.09% Sodium Azide as a preservative.

Purification: Processed to prevent minimal haemolysis and filtered to 0.2 μm.

Applications: Blocking sera.

Control sera.

Serum additive in dilution buffers.

Other applications not tested. Optimal dilutions are dependent on conditions and

should be determined by the user.

Add. Information: Source: Adult New Zealand Rabbits.

Storage: Store at 2-8°C for one month or in (aliquots) at -20°C for longer.

Avoid repeated freezing and thawing. Shelf life: one year from despatch.