

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES Phone: +1-888-267-4436 Fax: +1-301-340-8606

techsupport@origene.com

## OriGene Technologies GmbH

Schillerstr. 5 32052 Herford GERMANY Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info-de@origene.com

## AP09218BT-N Polyclonal Antibody to GFP (Ads. to Hu, Ms, Rt Serum Proteins) -

**Biotin** 

Alternate names: GFP-Tag, Green fluorescent protein

Quantity: 0.1 mg

Concentration: 1.0 mg/ml (by UV absorbance at 280 nm)

Background: Green fluorescence protein (GFP) is a 27 kDa protein derived from the jellyfish

Aequorea victoria, which emits green light (emission peak at a wavelenth of 509 nm)

when excited by blue light (excitation peak at a wavelenth of 395 nm).

Green Fluorescent Protein (GFP) has become an invaluable tool in cell biology research, since its intrinsic fluorescence can be visualized in living cells. GFP fluorescence is stable under fixation conditions and suitable for a variety of

applications.

GFP has been widely used as a reporter for gene expression, enabling researchers to visualize and localize GFP-tagged proteins within living cells without the need for

chemical staining.

Other applications of GFP include assessment of protein protein interactions through the yeast two hybrid system and measurement of distance between proteins through

fluorescence energy transfer (FRET) protocols.

GFP technnology has considerably contributed to a greater understanding of cellular

physiology.

Uniprot ID: <u>P42212</u> NCBI: <u>6100</u>

Host / Isotype: Rabbit / IgG

Immunogen: GST-Green Fluorescent Protein (GFP) fusion protein corresponding to the full length

amino acid sequence (246aa) derived from the jellyfish Aequorea victoria

**Format:** State: Lyophilized purified Ig fraction

**Purification:** Immunoaffinity Chromatography using Green Fluorescent Protein (Aequorea victoria) coupled to agarose beads followed by solid phase adsorption(s)

to remove any unwanted reactivities.

**Buffer System:** 0.02M Potassium Phosphate, 0.15M Sodium Chloride, pH 7.2, 10 mg/ml BSA (IgG and Protease free) as stabilizer and 0.01% (w/v) Sodium Azide as

preservative Label: Biotin

Molar Ratio: 10-20 BAC molecules per Rabbit IgG molecule

**Reconstitution:** Resuspend product with 0.1 ml of deionized water (or equivalent).

Applications: Suitable for Immunoblotting, ELISA, Immunohistochemistry, Immunomicroscopy as

well as other antibody based assays using streptavidin or avidin conjugates requiring

lot-to-lot consistency.
Recommended Dilutions:
ELISA: 1/10,000-1/50,000.

Western blot: 1/2,000-1/10,000.



**Biotin** 

Immunohistochemistry: 1/1,000-1/5,000.

Other applications not tested. Optimal dilutions are dependent on conditions and

should be determined by the user.

**Specificity:** wt, rGFP, eGFP.

Assay by Immunoelectrophoresis resulted in a single precipitin arc against anti-Rabbit Serum, anti-Biotin and purified and partially purified Green Fluorescent Protein

(Aequorea victoria).

No reaction was observed against Human, Mouse and Rat Serum Proteins.

Storage: Store vial at 2-8°C prior to restoration. Following restoration product can be stored

undiluted at 2-8°C for one month or (in aliquots) at -20°C or below.

Avoid repeated freezing and thawing.

Centrifuge product if not completely clear after standing at room temperature.

Shelf life: One year from despatch.

General Readings: 1. Bayer & Wilchek Methods in Enzymology 184; 138-160, 1990 (Conjugation).