

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES Phone: +1-888-267-4436 Fax: +1-301-340-8606

techsupport@origene.com

OriGene Technologies GmbH

Schillerstr. 5 32052 Herford GERMANY Phone: +49-5221-346

Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info-de@origene.com

AP09240CP-N Btb1 control peptide

Quantity: 50 μg

Concentration: 1.0 mg/ml (by weight)

Format: State: Liquid (sterile filtered) peptide

Purification: Greater than 95% specific peptide

Buffer System: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Applications: Control peptide should be used at 1.0 μg per 1.0 μl of antiserum per assay.

Other applications not tested. Optimal dilutions are dependent on conditions and

should be determined by the user.

Specificity: Control peptide for antibody AP09240PU-N only.

Storage: Store the antibody at -20°C.

Avoid repeated freezing and thawing. Shelf life: 6 month from despatch.