

OriGene Technologies Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES

Phone: +1-888-267-4436 Fax: +1-301-340-8606 techsupport@origene.com AP10264PU-N

OriGene EU

Acris Antibodies GmbH

Schillerstr. 5 32052 Herford GERMANY

Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info@acris-antibodies.com

Polyclonal Antibody to PGCC (C-term) - Aff - Purified

Alternate names: PGC-C, Particulate Guanylyl Cyclase C

Catalog No.: AP10264PU-N

Quantity: 0.1 mg

Concentration: 0.75 - 1.2 mg/ml (lot specific)

Background: Cyclic GMP (cGMP), a key messenger in several signal transduction pathways, the

intracellular levels of cGMP ate maintained by the activity of opposing enzymes: synthesizing gualylyl cyclases (GC) and hydrolyzing phosphodiesterases (PDEs). The synthesizing enzymes (GCs) are found in two forms: cytosolic (soluble) and membranebound (particulate), while they share similar structural characteristics, they differ in their mechanisms of physiological regulations. Most importantly, sGC contains a heme group and binds NO that activates the enzyme, while particulate GC is stimulated by natriuretic peptides. Plasma membrane forms of guanylyl cyclase have been shown to function as natriuretic peptide receptors. In response to G-protein couples receptor stimulation, the cGMP can be produced from GTP by either cytoplasmic, soluble guanylate cyclase (sGC) are heterodimers (a & b polypeptide chains), that are stimulated by nitric oxide and carbon monoxide or by particulate membrane-bound guanylyl cyclases which are activated by a complex mechanism by natriuretic peptides. Particulate GC (PGCs) have 7 different isoforms, PGC-A through PGC-G and are expressed in most tissues in isoform specific manner (See Table 1). There is significant structural homology among various PGCs, there is a large N-terminal extracellular domain (ECD), a single TMD and a large intracellular domain with protein kinase activity (KLD), a C-terminal catalytic domain (CD) and in between is a dimmerization domain (DD). The particulate GC-D, also known as olfactory gunylate cyclas, is specifically expressed in a subpopulation of olfactory sensory neurons. The expression of GC-D RNA is restricted to a small, randomly dispersed population of

neurons that is within a single topographic zone in the olfactory neuroepithelium and resembles the pattern of the more diverse seven-transmembrane-domain odorant

receptors (1). Particulate GC-D resembles two other GC expressed in retina (P-GC-E and F), but diverges considerable form other GCs. These observations suggest that GC-D may function directly in odor recognition or in modulating the sensitivity of a subpopulation of

sensory neurons to specific odors.

Host: Rabbit

Immunogen: Synthetic C-terminal peptide for particulate guanylyl cyclase C

Format: State: Liquid Ig fraction

Purification: Affinity chromatography **Buffer System:** Stabilization buffer

Applications: Western blot: > 1:500.

TUV HORD CERT COME!



AP10264PU-N: Polyclonal Antibody to PGCC (C-term) - Aff - Purified

Other applications not tested. Optimal dilutions are dependent on conditions and should

be determined by the user.

Specificity: This antibody detects PGCC at C-term.

Species Reactivity: Tested: Human, mouse, rat

Storage: Store (in aliquots) at -20 °C. Avoid repeated freezing and thawing.

Shelf life: one year from despatch.

General Readings: 1. Wedel B. J and Garbers D. L., Trend Endocrinol. Met. 9, 213-219; 1998

2. Hidaka Y, Matsumoto Y, Shimonishi Y. FEBS Lett 2002 526; 58-62 3. Deguch T., Amano E., Nakeane M. J. Neurochem. 27, 1027-1034, 1976.

4. Teunissen C et. al,. Brain Res. 891, 206-212; 2001.

