

OriGene Technologies, Inc.


9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES Phone: +1-888-267-4436 Fax: +1-301-340-8606 techsupport@origene.com

OriGene Technologies GmbH

Schillerstr. 5 32052 Herford GERMANY Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info-de@origene.com

AP23453SU-N Polyclonal Antibody to Domoic Acid - Serum

| Quantity: | 0.1 ml | |
|---------------|---|---|
| Host: | Rabbit | |
| Immunogen: | Immunization with Domoic Acid coupled to immunization with Domoic Acid coupled to | |
| Format: | State: Lyophilized Serum Reconstitution: Restore in aqua bidest to initial volume. | |
| Applications: | ELISA: 1/2,000. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user. | |
| Specificity: | Domoic Acid. | |
| Storage: | Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing. Shelf life: one year from despatch. | |
| Pictures: | Domoic Acid | 0 |



MS/20170526

1/1

For research and in vitro use only. Not for diagnostic or therapeutic work. Material Safety Datasheets are available at www.acris-antibodies.com or on request.