

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES Phone: +1-888-267-4436 Fax: +1-301-340-8606 techsupport@origene.com

# OriGene Technologies GmbH

Schillerstr. 5 32052 Herford **GERMANY** Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info-de@origene.com

#### Polyclonal Antibody to Porcine IgG (Fc specific) - Serum AP33086SU-N

Pig IgG, Pig Immunoglobulin G Alternate names:

Quantity: 1 ml

**Concentration:** Total protein and IgG concentration in the antiserum are comparable to those of

pooled rabbit serum. No foreign proteins added.

Host:

Highly purified normal IgG isolated from pooled swine serum. Immunogen:

Freund's complete adjuvant is used in the first step of the immunization procedure.

Format: State: Delipidated, heat inactivated lyophilized stable whole serum

**Preservatives:** None, as it may interfere with the antibody activity.

Stabilizers: None

Reconstitution: Restore with 1 ml sterile distilled water.

**Applications:** Immunoprecipitation.

> Can be used in precipitating techniques as immunoelectrophoresis and radial immunodiffusion to identify the presence of IgG in swine serum or other body fluids or to determine its concentration. To prepare an immunoadsorbent for the purification of swine IgG from serum or plasma. The antiserum is not intended for use in nonprecipitating antibody-binding or other highly sensitive techniques. For such uses labelled and unlabelled cytochemical grade reagents meeting the necessary additional specificity and performance requirements are available.

**Recommended Working Dilutions:** 

Immunoelectrophoresis: Use 2 μl or equivalent against 120 μl antiserum.

Double Radial Immunodiffusion (Ouchterlony): Use a rosette arrangement with 10 µl antiserum in a 3 mm diameter centre well and 2 µl serum samples (neat and diluted)

in 2 mm diameter peripheral wells.

Antibody titre: 1/64 when tested against pooled normal sheep serum in agar-block

immunodiffusion titration.

Other applications not tested. Optimal dilutions are dependent on conditions and

should be determined by the user.



## **Specificity:**

The reactivity of the antiserum is restricted to the Fc part of the IgG molecule. In immunoelectrophoresis and radial immunodiffusion, using various antiserum concentrations against normal swine serum a single precipitin line is obtained which shows a reaction of identity with the precipitin line obtained with purified IgG. No precipitation reaction is obtained with purified IgA, IgM, and IgG/Fab fragments.

Cross-reactivity: The antiserum does not cross-react with any other component of the swine immunoglobulin system. Inter-species cross-reactivity is a normal feature of antibodies to immunoglobulins, since Ig of different species frequently share antigenic determinants. Cross-reactivity of this antiserum has not been tested in detail, however in immunoelectrophoresis a cross-reaction with serum of bovine, goat and human has been observed.

Species: Pig.

Other species not tested.

## Add. Information:

**Adsorption:** Immunoaffinity adsorbed using insolubilized antigens as required, to eliminate antibodies crossreacting with other components of the immunoglobulin system or reacting with other serum proteins. Special attention is given to the removal of antibodies to common Ig/Fab. The use of insolubilized adsorption antigens prevents the presence of excess adsorbent protein or immune complexes in

the antiserum.

### Storage:

Store lyophilized at 2-8°C for 6 months or at -20°C long term.

After reconstitution store the antibody undiluted at 2-8°C for one month

or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing. Shelf life: one year from despatch.