

OriGene Technologies Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES

Phone: +1-858-888-7900 Fax: +1-858-888-7904 <u>US-info@acris-antibodies.com</u> AR00140PU-N

OriGene EU

Acris Antibodies GmbH

Schillerstr. 5 32052 Herford GERMANY

Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info@acris-antibodies.com

S. aureus Teichoic Acid, Staphylococcus aureus (Wood Strain) - Purified

Catalog No.: AR00140PU-N

Quantity: 1 ml

Concentration: 2,0 mg/ml (RID) Titer: 1:16 endpoint dilution (RID)

Species: S. aureus

Source: Staphyloccus aureus (Wood Strain)

Format: State: Liquid purified Ig

Purity: >95% pure (SDS-PAGE) Buffer System: PBS, pH 7.2

Description: Teichoic Acid from Staphyloccus aureus (Wood Strain)

Inactivation: Sterile filtered. Sterility test: No growth

Storage: Store the protein undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

Avoid repeated freezing and thawing. Shelf life: one year from despatch.

Caution: All materials should be handled as if potentially infectious. Generally accepted laboratory

practices appropriate for infectious materials should be employed when handling this

product.

