

## OriGene Technologies Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES

Phone: +1-858-888-7900 Fax: +1-858-888-7904 US-info@acris-antibodies.com

## AR10521PU-L OriGene EU

## **Acris Antibodies GmbH**

Schillerstr. 5 32052 Herford GERMANY

Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info@acris-antibodies.com

## Dengue Virus Polyvalent type 1-4 - Purified

Catalog No.: AR10521PU-L

Quantity: 1 mg

Concentration: 0.79 mg/ml

Background: Caused by one of four closely related virus serotypes of the genus Flavivirus, family

Flaviviridae, each serotype is sufficiently different that there is no cross-protection and epidemics caused by multiple serotypes (hyperendemicity) can occur. In cell culture experiments and mice Morpholino antisense oligos have shown specific activity against

Dengue virus.

Source: E. coli

Format: Purity: >95% pure as determined by 10% PAGE (coomassie staining).

Purification Method: Purified by proprietary chromatographic technique.

Buffer System: PBS, pH 7.4

**Description:** Polyvalent Dengue antigen is a 22kDa protein that contains 4 Dengue antigen subtypes,

produced in *E.coli* and fused to a 6xHis Tag.

The immunoreactive peptide span from the viral envelope to the non-functional region. The polyvalent Dengue antigen consists of 25% recombinant antigen from each of the Dengue

subtypes in their composition.

The Polyvalent Dengue antigen detects all four dengue subtypes.

The sensitivity and specificity could reach approximately 93% and 96% in Dengue IgG and

IgM ELISA assay.

Storage: Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

Avoid repeated freezing and thawing. Shelf life: one year from despatch.