

OriGene Technologies Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES Phone: +1-858-888-7900 Fax: +1-858-888-7904 <u>US-info@acris-antibodies.com</u>

AR10696PU-L OriGene EU

Acris Antibodies GmbH Schillerstr. 5 32052 Herford GERMANY Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info@acris-antibodies.com

VZV / HHV-3 Glyco E (48-135) - Purified

Catalog No.:	AR10696PU-L
Quantity:	1 mg
Concentration:	1.0 mg/ml
Background:	VZV is closely related to the herpes simplex viruses(HSV), sharing much genomehomology. The known envelope glycoproteins (gB, gC, gE, gH, gI, gK, gL) correspond with those in HSV, however there is no equivalent of HSV gD. VZV virons are spherical and 150-200 nm in diameter. Their lipidenvelope encloses the nucleocapsid of 162 capsomeres arranged in a hexagonalform. Its DNAis a single, linear, double-stranded molecule, 125,000 nt long. The virus is very susceptible to disinfectants, notably sodium hypochlorite. Within the body it can be treated by a number of drugs and therapeutic agents including aciclovir, zoster- immune globulin(ZIG), and vidarabine.
Source:	E. coli
Format:	Purity: >95.0% pure as determined by 10.0% PAGE (coomassie staining). Purification Method: Sepharose-Derived Purification. Buffer System: 25mM Tris-HCl, 1mM EDTA, 50% Glycerol
Applications:	Antigen in ELISA and Western blots, excellent antigen for detection of VZV with minimal specificity problems. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Description:	<i>E.coli</i> derived recombinant protein contains the VZV gE immunodominant regions. Specificity: Immunoreactive with sera of VZV-infected individuals.
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. Shelf life: one year from despatch.



MS/20120905