

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES Phone: +1-888-267-4436 Fax: +1-301-340-8606 techsupport@origene.com

OriGene Technologies GmbH

Schillerstr. 5 32052 Herford GERMANY Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info-de@origene.com


AR50536PU-S Human MOCS2 (1-188, His-tag) - Purified

Alternate names:	MCBPE, MOCO1, MOCS2A, MOCS2B, MPTS, molybdenum cofactor synthesis 2
Quantity:	20 μg
Concentration:	0.25 mg/ml (determined by Bradford assay)
Background:	MOCS2 is a heterotetrameric synthase composed of two small (MOCS2A) and two large (MOCS2B) subunits. The large and small subunits of molybdopterin synthase are both encoded from this gene by overlapping open reading frames. MOCS2 functions in the second step of the synthesis of molybdenum cofactor or molybdopterin (MPT). It catalyzes the formation of MPT from precursor Z by incorporating a dithiolene functional group.
Uniprot ID:	<u>096007</u>
NCBI:	<u>NP_004522</u>
GenelD:	<u>4338</u>
Species:	Human
Source:	E. coli
Format:	State: Liquid purified protein Purity: >90% by SDS - PAGE Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 50% glycerol,1mM DTT
Description:	Recombinant human MOCS2 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques. AA Sequence: MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSMSSL EISSSCFSLE TKLPLSPPLV EDSAFEPSRK DMDEVEEKSK DVINFTAEKL SVDEVSQLVI SPLCGAISLF VGTTRNNFEG KKVISLEYEA YLPMAENEVR KICSDIRQKW PVKHIAVFHR LGLVPVSEAS IIIAVSSAHR AASLEAVSYA IDTLKAKVPI WKKEIYEESS TWKGNKECFW ASNS Molecular weight: 25 kDa (224aa) confirmed by MALDI-TOF
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing. Shelf life: one year from despatch.
General Readings:	Leimkuhler S., et al. (2005) Hum Genet. 117: 565-570Reiss J., et al. (1999) Hum Genet . 64: 706-711.

For research and in vitro use only. Not for diagnostic or therapeutic work. Material Safety Datasheets are available at www.acris-antibodies.com or on request.

OG/20121030 1/2

Pictures:



15% SDS-PAGE (3ug)

For research and in vitro use only. Not for diagnostic or therapeutic work. Material Safety Datasheets are available at www.acris-antibodies.com or on request.