

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850 UNITED STATES Phone: +1-888-267-4436 Fax: +1-301-340-8606 techsupport@origene.com

OriGene Technologies GmbH

Schillerstr. 5 32052 Herford GERMANY Phone: +49-5221-34606-0 Fax: +49-5221-34606-11 info-de@origene.com

BA1089 Human HDL 2 Subfraction - Purified

Alternate names: Quantity: Concentration: Species:	HDL, Lipoprotein High Density 2 Subfraction 1 mg 1.0 mg/ml (BCA with SDS)(lot specific) Human
Source:	Plasma, (Human Normolipidemic Plasma)
Format:	 State: Liquid purified fraction. Purity: Sequential Density Ultracentrifugation at density 1.063-1.120. Devoid of VLDL, IDL, LDL, Lp[a] and HDL3. Buffer System: 100mM Sodium Chloride, 10mM Tris, 1mM EDTA, pH 7.4 Preservatives: 1mM Sodium Azide
Description:	Human High Density Lipoprotein-2 (HDL2) Subfraction.
Storage:	Store at 2-8°C. DO NOT FREEZE!! Shelf life: one year from despatch.
Caution:	All human source materials have tested negative for HIV1, HIV2, HCV antibodies and HBsAg. No test guarantees a product to be non-infectious. Therefore, all material derived from human fluids or tissues should be considered as potentially infectious.

For research and in vitro use only. Not for diagnostic or therapeutic work. Material Safety Datasheets are available at www.acris-antibodies.com or on request.

1/1