



OriGene Technologies Inc.

9620 Medical Center Drive, Ste 200  
Rockville, MD 20850  
UNITED STATES  
Phone: +1-858-888-7900  
Fax: +1-858-888-7904  
[US-info@acris-antibodies.com](mailto:US-info@acris-antibodies.com)

BP046

OriGene EU

Acris Antibodies GmbH

Schillerstr. 5  
32052 Herford  
GERMANY  
Phone: +49-5221-34606-0  
Fax: +49-5221-34606-11  
[info@acris-antibodies.com](mailto:info@acris-antibodies.com)

## Polyclonal Antibody to Amphibian Bombesin

|                       |   |
|-----------------------|---|
| <b>Catalog No.:</b>   | BP046   |
| <b>Quantity:</b>      | 0.25 ml   |
| <b>Concentration:</b> | conc. lot specific  |
| <b>Host:</b>          | Rabbit  |
| <b>Immunogen:</b>     | BSA conjugate   |
| <b>Applications:</b>  | Immunohistochemistry on paraffin embedded sections: 1/400-1/2000 (ABC).<br>Immunofluorescence: 1/30-1/300. Positive control: Human foetal lung. Fixatives: Bouins is recommended. Other applications not tested. Optimal dilutions of this antibody are dependent on conditions and should be determined by the user.<br>Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user. |
| <b>Specificity:</b>   | 1 nmol synthetic amphibian bombesin per ml diluted antiserum abolishes staining.  |
| <b>Storage:</b>       | Aliquot and store at -20°C. Avoid repeated freezing and thawing. Shelf life: on   |