



OriGene Technologies Inc.

9620 Medical Center Drive, Ste 200
Rockville, MD 20850
UNITED STATES
Phone: +1-858-888-7900
Fax: +1-858-888-7904
US-info@acris-antibodies.com

BP1012

OriGene EU

Acris Antibodies GmbH
Schillerstr. 5
32052 Herford
GERMANY
Phone: +49-5221-34606-0
Fax: +49-5221-34606-11
info@acris-antibodies.com

Polyclonal Antibody to Clostridium Botulinus (Toxin)

Catalog No.:	BP1012
Quantity:	1 ml
Concentration:	4.22 mg/ml (OD280, lot specific)
Host:	Rabbit
Immunogen:	Enzymatically denatured Clostridium botulinus toxin (Subtype A).
Applications:	ELISA (Titre is dependant on the quality of the antigen. We recommended that this antibody is titred in the end-users system using a dilution range of 1/1000 - 1/100K. In a simple ELISA with 1 µg/ml non-denatured Botulinus toxin in the coating buffer, the end point was observed at an antibody dilution of 1/40K). Other applications not tested. Optimal dilutions of this antibody are dependent on conditions and should be determined by the user. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	Clostridium Botulinus (Toxin)
Storage:	Aliquote and store the antibody at -20°C. Avoid repeated freezing and thawing. Shelf life: on