



OriGene Technologies Inc.

9620 Medical Center Drive, Ste 200  
Rockville, MD 20850  
UNITED STATES  
Phone: +1-858-888-7900  
Fax: +1-858-888-7904  
[US-info@acris-antibodies.com](mailto:US-info@acris-antibodies.com)

BP1014

OriGene EU

Acris Antibodies GmbH

Schillerstr. 5  
32052 Herford  
GERMANY  
Phone: +49-5221-34606-0  
Fax: +49-5221-34606-11  
[info@acris-antibodies.com](mailto:info@acris-antibodies.com)

## Polyclonal Antibody to Clostridium perfringens (Enterotoxin A)

<b>Catalog No.:</b>	BP1014
<b>Quantity:</b>	1 ml
<b>Concentration:</b>	7.7 mg/ml
<b>Host:</b>	Goat
<b>Immunogen:</b>	native Clostridium perfringens (enterotoxin A)
<b>Applications:</b>	ELISA (to be determined; In a simple ELISA with 1 µg/ml antigen coated, the end point dilution was observed at an antibody dilution of >1/100,000). RID/IEP (1/2 - 1/5). Western blotting (1/2 - 1/5). IHC (1/20 - 1/200; positive binding was observed in human stools known to contain enterotoxin from C1.). Other applications not tested. Optimal dilutions of this antibody are dependent on conditions and should be determined by the user. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
<b>Specificity:</b>	Clostridium perfringens (enterotoxin A). Cross reactivity: None known.
<b>Storage:</b>	Aliquot and store the antibody at -20°C. Avoid repeated freezing and thawing. Shelf life: on