



OriGene Technologies Inc.

9620 Medical Center Drive, Ste 200
Rockville, MD 20850
UNITED STATES
Phone: +1-858-888-7900
Fax: +1-858-888-7904
US-info@acris-antibodies.com

BP400

OriGene EU

Acris Antibodies GmbH
Schillerstr. 5
32052 Herford
GERMANY
Phone: +49-5221-34606-0
Fax: +49-5221-34606-11
info@acris-antibodies.com

Polyclonal Antibody to Crystallin beta H - Purified

Alternate names:	beta H Crystallin
Catalog No.:	BP400
Quantity:	1 ml
Concentration:	17.4 mg/ml
Host:	Sheep
Immunogen:	Native Crystallin, from eye lens
Format:	State: Liquid Ig fraction Purification: Prepared by Ammonium Sulphate fractionation Buffer System: PBS, pH 7.2, 0.09 % Sodium Azide
Applications:	ELISA. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	This antibody detects Beta-H-Crystallin 100 %, beta-L-Crystallin 40 %, alpha-Crystallin 10 %. Species: Bovine. Other species not tested.
Storage:	Store the antibody at -20 °C. Can be shipped at 2 - 8 °C. Avoid repeated freezing and thawing. Shelf life: One year from despatch.