



OriGene Technologies Inc.

9620 Medical Center Drive, Ste 200
Rockville, MD 20850
UNITED STATES
Phone: +1-888-267-4436
Fax: +1-301-340-8606
techsupport@origene.com

BP646

OriGene EU

Acris Antibodies GmbH
Schillerstr. 5
32052 Herford
GERMANY
Phone: +49-5221-34606-0
Fax: +49-5221-34606-11
info@acris-antibodies.com

Polyclonal Antibody to Neomycin

Catalog No.:	BP646
Quantity:	0.5 ml
Concentration:	29.5 mg/ml (OD280, lot specific)
Host:	Rabbit
Immunogen:	Neomycin KLH-Konjugat
Applications:	ELISA (1/30 for hap-BSA, 50%B; 1/20,000 for hap-BSA, EPT): Use at 1/100 to 1/1000 final dilution in simple ELISA with overnight incubation. Under above conditions can be used to specifically detect Neomycin in the range 0.1 to 100pM. This antibody may also be used to measure Tobramycin (wide assay range) with altered antibody dilutions. Other applications not tested. Optimal dilutions of this antibody are dependent on conditions and should be determined by the user. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	Neomycin 100%, Tobramycin 0.07%.
Storage:	Store the antibody at 4-8°C for one month or at -20°C for longer. Avoid repeated freezing and thawing. Shelf life: on